



| | |
|--------------------------------------|-----------------------|
| Chromatin Immunoprecipitation (ChIP) | |
| Document #: LIBPR.0015 | Supersedes: Version 7 |
| Version: 8 | Page 1 of 17 |

Non Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

Chromatin Immunoprecipitation (ChIP)

I. Purpose

To perform Chromatin Immunoprecipitation (ChIP) on crosslinked cell lines and tissues.

II. Scope

All procedures are applicable to the BCGSC Library Core group and Library Technology Development group.

III. Policy

This procedure will be controlled under the policies of the Genome Sciences Centre, as outlined in the Genome Sciences Centre High Throughput Production Quality Manual (QM.0001). Do not copy or alter this document. To obtain a copy see a QA associate.

IV. Responsibility

It is the responsibility of all personnel performing this procedure to follow the current protocol. It is the responsibility of the Group Leader to ensure personnel are trained in all aspects of this protocol. It is the responsibility of Quality Assurance Management to audit this procedure for compliance and maintain control of this procedure.

V. References

| Document Title | Document Number |
|----------------|-----------------|
| N/A | N/A |

VI. Related Documents

| Document Title | Document Number |
|---|-----------------|
| Operation and Maintenance of the Agilent 2100 Bioanalyzer for DNA samples | LIBPR.0017 |
| Quantifying DNA Samples using the Qubit Fluorometer | LIBPR.0030 |
| qPCR of ChIP Samples | LIBPR.0075 |

VII. Safety

All Laboratory Safety procedures will be complied with during this procedure. The required personal protective equipment includes a laboratory coat and gloves. See the material safety data

| Chromatin Immunoprecipitation (ChIP) | |
|--------------------------------------|-----------------------|
| Document #: LIBPR.0015 | Supersedes: Version 7 |
| Version: 8 | Page 2 of 17 |

Non Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

sheet (MSDS) for additional information. **Steps that involve the use of PCI should be performed in the fume hood.**

VIII. Materials and Equipment

| Name | Supplier | Number: # | Model or Catalogue # |
|---|-----------------------|-----------------------|----------------------|
| Small size safetouch nitrile gloves | Ultident | 296359683 | ✓ |
| wet ice | In house | N/A | N/A |
| Ice bucket | Fisher | 11-675-58 | ✓ |
| 1.5 mL Microtubes | Diamed | PRE150-B | ✓ |
| 15mL Conical Tubes | BD Falcon | 352097 | ✓ |
| 50mL Conical Tubes | BD Falcon | 352070 | |
| Gilson P10 pipetman | Mandel | GF-44802 | ✓ |
| Gilson P20 pipetman | Mandel | GF23600 | ✓ |
| Gilson P200 pipetman | Mandel | GF-23601 | ✓ |
| Gilson P1000 pipetman | Mandel | GF-23602 | ✓ |
| Diamond Filter tips DFL10 (10 Tipacks of 96 racked filter tips) | Mandel Scientific | GF-F171203 | ✓ |
| Diamond Filter tips DFL30 (10 Tipacks of 96 racked filter tips) | Mandel Scientific | GF-F171303 | ✓ |
| Diamond Filter tips DFL200 (10 Tipacks of 96 racked filter tips) | Mandel Scientific | GF-F171503 | ✓ |
| Diamond Filter tips DFL1000 (10 Tipacks of 96 racked filter tips) | Mandel Scientific | GF-F171703 | ✓ |
| Galaxy mini-centrifuge | VWR | 37000-700 | ✓ |
| Large Kimwipes (Kimberly Clark/Kimtech) | Fisher Scientific | 06-666-1A | ✓ |
| Black ink permanent marker pen | VWR | 52877-310 | ✓ |
| Small Autoclave waste bags 10"X15" | Fisher Scientific | 01-826-4 | ✓ |
| 10ml serological pipettes | Costar | CS004488 | ✓ |
| Portable Pipet Aid, Multispeed XP, rechargeable | Fisher Scientific | 13-681-15E | ✓ |
| Ultra Pure Water (Rnase/Dnase free) | Invitrogen | 10977-023 | ✓ |
| Anhydrous Ethyl Alcohol (100% Ethanol) | Commercial Alcohols | People Soft ID: 23878 | ✓ |
| DNA away | Molecular Bioproducts | 7010 | ✓ |
| nProtein A Sepharose 4 Fast Flow | GE Healthcare | 17-5280-01 | ✓ |
| nProtein G Sepharose 4 Fast Flow | GE Healthcare | 17-0618-01 | ✓ |
| 1M Tris Solution, pH 8.0 | Ambion | AM9856 | ✓ |
| Triton® X-100, laboratory grade | Sigma | X100-100ML | ✓ |
| Deoxycholic acid, sodium salt | Fisher Scientific | AC21859-0250 | ✓ |
| 20% SDS Solution | Ambion | AM9820 | ✓ |
| 5M Sodium Chloride | Ambion | AM9760G | ✓ |
| 0.5M EDTA, pH8.0 | Ambion | AM9260G | ✓ |
| Complete-Mini EDTA-free Protease Inhibitor Cocktail Tablet EasyPack | Roche Diagnostics | 04 695 159 001 | ✓ |



| Chromatin Immunoprecipitation (ChIP) | |
|--------------------------------------|-----------------------|
| Document #: LIBPR.0015 | Supersedes: Version 7 |
| Version: 8 | Page 3 of 17 |

Non Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

| | | | | |
|--|-------------------|----------------|-----|-----|
| Complete-Mini Protease Inhibitor Cocktail Tablet EasyPack | Roche Diagnostics | 04 693 124 001 | | ✓ |
| Rnase-free Non-stick tubes, 1.5mL, | Ambion | AM12450 | | ✓ |
| Centrifuge, Eppendorf 5417R, refrigerated high-speed, 115V | Fisher Scientific | 5417 R | ✓ | |
| Deoxyribonucleic acid, single stranded from salmon testes | Sigma | D7656-1ML | | ✓ |
| Bovine Serum Albumin (BSA), Buffers, 10mg/ml | NEB | B9001S | | ✓ |
| Non-stick 2ml tubes | Ambion | AM12475 | | ✓ |
| Syringe, 1ml, Gauge: 26 | Fisher Scientific | 333-309625 | | ✓ |
| Syringe, 3ml, Gauge: 21 | Fisher Scientific | 14-826-84 | | ✓ |
| LabQuake Shaker/Rotator with Clips | Barnstead | 415110 | | ✓ |
| Sonic Dismembrator 550 (cup horn) | Fisher Scientific | Discontinued | N/A | N/A |
| Foam tube holder | Ambion Sample | N/A | N/A | N/A |
| QIAquick PCR Purification Kit, 50rxn | Qiagen | 28104 | | ✓ |
| Proteinase K (20µg/µL) | Ambion | AM2546 | | ✓ |
| 2mL phase lock gel heavy tubes | Inter Medico | 5P2302830 | | ✓ |
| Phenol/Chloroform/Isoamyl Alcohol | Fisher Scientific | BP1752-400 | | ✓ |
| Sodium Acetate, 3M, pH 5.5 | Ambion | AM9740 | | ✓ |
| | Mandel Scientific | GF-F164001 | | ✓ |
| DistriMan Pipettor | | | | ✓ |
| 96 Place Microfuge Tube Rack | WVR | CBGTR-096 | | ✓ |
| 0.5mL Non-stick Microfuge Tube | Ambion | AM12350 | N/A | N/A |
| VX-100 Vortex Mixer | Rose Scientific | S-0100 | ✓ | |
| Sodium bicarbonate | Sigma | S7277-250G | | ✓ |
| Rnase, Dnase-free, 500µg | Roche | 11119915001 | | ✓ |
| Eppendorf Thermomixer 1.5 mL | Eppendorf | 21516-166 | ✓ | |
| Parafilm | Fisher Scientific | 13-374-12 | | ✓ |
| Mussel Glycogen | Roche | 10 901 393 001 | | ✓ |
| Spin-X Filter Tube | Fisher Scientific | CS008160 | | ✓ |
| | Mandel Scientific | | | ✓ |
| DistriMan Tips (12.5mL capacity) | | GF-F164150 | | ✓ |
| Liquid Nitrogen | Praxair | NI M-FILLTR | | ✓ |
| H3K4me1 | Diagenode | C15410037 | | ✓ |
| H3K4me3(C42D8) Rb mAb | Cell Signaling | 9751S | | ✓ |
| H3K9me3 | Diagenode | C15410056 | | ✓ |
| H3K27me3 | Diagenode | C15410195 | | ✓ |
| H3K36me3 | Abcam | ab9050 | | ✓ |
| H3K27ac | Abcam | ab4729 | | ✓ |
| Closed System Tissue Grinder | Fisher Scientific | 02-542-09 | | ✓ |
| 1M DTT | Invitrogen | P2325 | | ✓ |
| Heated Lids | In House | n/a | | |

IX. Procedure

6th floor ChIP Room



| | |
|--------------------------------------|-----------------------|
| Chromatin Immunoprecipitation (ChIP) | |
| Document #: LIBPR.0015 | Supersedes: Version 7 |
| Version: 8 | Page 4 of 17 |

Non Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

1. Retrieval of reagents, sample(s) and equipment preparation (Day 1)

***** Consult with supervisor on the status of the sample. If working with cells it is important to know the cell species, amount of cells, vessel the cells are in, if the cells are crosslinked or not, and how long the cells were crosslinked for if they already are. If working with tissue, it is important to know the tissue species, approximate weight of tissue, if the tissue has been crosslinked already and if so, how long the tissue was crosslinked for. Also regardless of the sample type, it is important to clarify as to approximately how many IPs are to be set up and with which antibodies. This helps determine how many beads to block.*****

- 1.1. Put on a clean pair of gloves and lab coat.
- 1.2. Wipe down the work bench, small equipment, and ice bucket with DNA Away, dH₂O and 70% Ethanol.
- 1.3. Change gloves.
- 1.4. Retrieve fresh ice and all required reagents.
- 1.5. Print out the ChIP worksheets located at R:\Library Core\Work Sheets and Calculators\Epigenomics\ChIP

2. Lysing of cells or tissue

- 2.1. Dissolve 1 Complete-Mini PIC tablet in 10mL of ChIP Lysis Buffer and mix by inverting. Keep the mix at room temperature.
- 2.2. If the sample(s) is a cell line or primary cells, proceed to Step 2.3. If the sample(s) is tissue, proceed to Step 2.4.
- 2.3. Based on the cell count, add the following amounts of ChIP Lysis Buffer + PIC to the cell pellet to resuspend. Do not invert the tube when resuspending the pellet. Proceed to Step 2.5.

| Cell Count | ChIP Lysis Buffer + PIC to add |
|--|--------------------------------|
| <1x10 ⁶ cells | 500µL |
| 1x10 ⁶ - 10x10 ⁷ cells | 1mL |
| ≥10x10 ⁷ cells | 2mL |



| | |
|--------------------------------------|-----------------------|
| Chromatin Immunoprecipitation (ChIP) | |
| Document #: LIBPR.0015 | Supersedes: Version 7 |
| Version: 8 | Page 5 of 17 |

Non Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

- 2.4. For tissue that was crosslinked by the collaborator, it is important to physically breakdown the tissue a bit before and during lysis. After adding lysis buffer (1-2mL) to the tissue, transfer the sample to the closed system tissue grinder tube. Consult with supervisor if unsure what volume of lysis buffer to add. The tissue sample should appear as a 'slurry' in lysis buffer. Securely close the lid. Ensure not to touch the tip of the grinder to prevent any contamination of the sample. Twist the handle of the grinder, to grind and break up the tissue along the sides of the tube. DO NOT plunge the handle up and down. The lysis buffer contains SDS and foam will generate. Once the tissue looks fully broken down, open the system and aspirate the sample out using a wide bore P1000 tip. Consult with supervisor if unclear. Proceed to Step 2.5.
- 2.5. Incubate on ice for 30 minutes. Beads can be prepped during this incubation (See Step 3).
- 2.6. Transfer a maximum of 1mL of lysed sample into a 2mL non-stick conical tube. Use more tubes if there is more than 1mL of sample.
- 2.7. Pass the sample through a chilled 1mL syringe (with 26 gauge) 5-6 times. If the sample is tissue, it is sometimes difficult to pass the lysed material through a 26 gauge needle. Pass through a 21 gauge needle.
- 2.8. Spin at 5000rpm for 10 minutes at 4°C.
- 2.9. Carefully and slowly remove supernatant as the pellet can be aspirated quite easily. The pellet will appear almost jelly-like.
- 2.10. Resuspend each pellet with ChIP Lysis Buffer + PIC to a final volume of 250uL per pellet. If the chromatin pellet is very large, the volume can be increased up to 300uL with ChIP Lysis Buffer + PIC.
- 2.11. Resuspend the pellet by finger vortexing and pulse spinning.
- 2.12. Proceed to Step 4.

3. Preparation of Beads

- 3.1. Calculate the amount of beads (1:1 mixture of Protein A Sepharose and Protein G Sepharose beads) required for IP set-up and preclearing the chromatin.

40µL of beads per sample (for preclearing) + (20µL of beads x number of IPs)

e.g. 2 samples and 16 IPs



| | |
|--------------------------------------|-----------------------|
| Chromatin Immunoprecipitation (ChIP) | |
| Document #: LIBPR.0015 | Supersedes: Version 7 |
| Version: 8 | Page 6 of 17 |

Non Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

$(40\mu\text{L} \times 2 \text{ samples}) + (20\mu\text{L} \text{ of beads} \times 16 \text{ IPs}) = 400\mu\text{L} \text{ of A/G Sepharose Bead Mix}$

Note: When pipetting out beads, use a p200 tip that has the end of the tip cut off for a wider opening. All subsequent handling of beads should be performed with p200 tips that have the end of the tip cut off.

- 3.2. Resuspend the stock of Protein G Sepharose Beads and Protein A Sepharose Beads by inverting the bottles several times.
- 3.3. Add Protein A Sepharose Beads and Protein G Sepharose beads together in a 1:1 ratio into a 1.5mL non-stick tube. See below using the example in Step 2.1.

$400\mu\text{L} \text{ of A/G Sepharose Bead Mix} = 200\mu\text{L} \text{ of Protein A Sepharose Beads} + 200\mu\text{L} \text{ of Protein G Sepharose Beads}$

- 3.4. Dissolve 1 tablet of Complete-Mini EDTA-free PIC into 10ml of IP Buffer and mix by inverting. Do not vortex the buffer. Keep the mix at 4°C or on ice.
- 3.5. Add 1mL of the IP Buffer and PIC mix to beads and wash by inverting several times.
- 3.6. Pellet the beads by centrifuging for 2 minutes at 4000rpm at 4°C.
- 3.7. Place beads on a rack and allow 30 seconds for the beads to fully settle.
- 3.8. Remove the supernatant with a 200µL pipette. Avoid disturbing the beads when discarding supernatant.
- 3.9. Repeat wash with 1ml of the IP Buffer and PIC mix.
- 3.10. Pellet beads by centrifuging for 2 minutes at 4000rpm at 4°C.
- 3.11. Place beads on a rack and allow 30 seconds for the beads to fully settle.
- 3.12. Remove the supernatant with a 200uL pipette. Avoid disturbing the beads when discarding supernatant.
- 3.13. Add 1mL IP Buffer and PIC to the beads and then add:

5µL of 10mg/mL Salmon Sperm DNA/100µL bead mix
13µL of 10mg/mL BSA/100µL bead mix



| | |
|--------------------------------------|-----------------------|
| Chromatin Immunoprecipitation (ChIP) | |
| Document #: LIBPR.0015 | Supersedes: Version 7 |
| Version: 8 | Page 7 of 17 |

Non Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

Adjust based on the same ratio if different volume of beads is used.

- 3.14. Parafilm the tube and rotate on rotator for 3hrs at 4°C.
- 3.15. Pellet beads as described above in Step 3.10-3.12.
- 3.16. Resuspend beads 1:1 with IP Buffer and PIC mix to make a 50% slurry.
- 3.17. Store slurry on ice until ready to use.

4. Sonication

Note: When performing sonication, wear protective ear phones.

- 4.1 Connect either end of the coaxial High Frequency Cable to the SHV connector on to the rear panel of generator and connect other end of cable to the SHV connector on the convertor. Push the connectors on and turn the chrome rings $\frac{1}{4}$ turn to secure the connectors.
- 4.2 Mount the convertor and the cup horn onto a laboratory stand with large clamps. Do not hold or clamp the convertor by the black portion or by the horn itself, only support the convertor by clamping around the broad (chrome-plated) section.
- 4.3 Put the two black screws into the cup horn.
- 4.4 Turn the generator on by pressing "I" on Power Switch. The switch will illuminate when the power is on.
- 4.5 Adjust the amplitude by turning the dial to "7".
- 4.6 Press CLEAR to select Programmed Mode.
- 4.7 Press PROG/DATA to select programmed mode. The display will show the program screen and memory location.
- 4.8 Press PROG key a second time to begin programming the selected memory location.
- 4.9 Enter the process time: 00:10:00 (10 minutes) and press enter.



| | |
|--------------------------------------|-----------------------|
| Chromatin Immunoprecipitation (ChIP) | |
| Document #: LIBPR.0015 | Supersedes: Version 7 |
| Version: 8 | Page 8 of 17 |

Non Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

Note: Sonication time needed might be different for different samples, optimization might be needed. Consult supervisor.

4.10 Enter pulse "ON" time: 00:30:00 (30 seconds) and press enter.

4.11 Enter pulse "OFF" time: 00:30:00 (30 seconds) and press enter.

4.12 Secure sample (2.0mL tube) in the cup horn by placing the tube in an adapter manufactured by the GSC Engineering Group. Ensure the bottom of the tube is directly over the hole on the radiating surface of the cup horn.

4.13 Fill the cup horn with water and a little bit of ice until the water level is level with the sample in the 2mL tube.

4.14 Press start to begin sonicating the sample.

4.15 Set a timer to count up. Sonication runs for 20 minutes.

4.16 When the sonicator is in the pulse "OFF" time, aspirate water out and replenish with more ice to maintain the water level as in Step 4.13.

4.17 Once the program finishes, it will stop on its own. The program timer will automatically reset itself.

4.18 Sonicate the remaining samples as described above.

4.19 Once sonication is done, turn the amplitude dial to "0". Press "O" on the Power Switch to turn the generator off.

4.20 Disconnect the cable from SHV connector on the convector.

4.21 Spin the samples at 14,000rpm for 12 min at 4°C.

4.22 Transfer the supernatant to a fresh 1.5mL non-stick tube labeled as sonicated chromatin. Pool the supernatant from the same crosslinked sample together.

4.23 Aliquot 5uL of chromatin into a new 1.5mL microfuge tube with 195uL of Ultrapure water for reverse-crosslinking QC.

4.24 Store the remainder of the sonicated chromatin at 4°C.



| | |
|--------------------------------------|-----------------------|
| Chromatin Immunoprecipitation (ChIP) | |
| Document #: LIBPR.0015 | Supersedes: Version 7 |
| Version: 8 | Page 9 of 17 |

Non Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

5 Reverse-Crosslinking for sonicated chromatin QC check

- 5.1 Reverse crosslink by adding 8 μ L of 5M NaCl and 10 μ L of 0.5 μ g/ μ L RNase (Roche).
- 5.2 Incubate at 68 $^{\circ}$ C for 30 minutes on the thermomixer.
- 5.3 Briefly spin down the tube and let cool to room temperature. Change the temperature setting on the thermomixer to 42 $^{\circ}$ C.
- 5.4 Add 1 μ L of Proteinase K (20 μ g/ μ L).
- 5.5 Incubate at 42 $^{\circ}$ C for 30 minutes on the thermomixer.
- 5.6 Spin down 2mL phase lock gel tubes at 14,000rpm for 1 min at room temperature.
- 5.7 Add sample (~200 μ L) to the 2mL phase lock gel tube and add 200 μ L PCI to the sample. Mix by inverting.
- 5.8 Centrifuge at 14,000rpm for 5 min at room temperature.
- 5.9 Transfer the supernatant to a fresh 1.5mL tube. Add 1 μ L of 10mg/mL mussel glycogen, 20 μ L of 3M NaOAc pH5.5, and 500 μ L of ice cold 100% EtOH. Vortex to mix completely and place at -20 $^{\circ}$ C for at least 20 minutes.
- 5.10 Centrifuge at 14,000rpm for at least 30 minutes at 4 $^{\circ}$ C.
- 5.11 Carefully remove and discard supernatant. Do not disturb pellet.
- 5.12 Wash pellet with 1mL 70% EtOH.
- 5.13 Spin at 14,000rpm for 5 min at 4 $^{\circ}$ C.
- 5.14 Remove and discard supernatant.
- 5.15 Repeat Steps 5.12 to Step 5.14.
- 5.16 Allow pellet to dry.
- 5.17 Resuspend pellet in 5 μ L of ultrapure water.



| | |
|--------------------------------------|-----------------------|
| Chromatin Immunoprecipitation (ChIP) | |
| Document #: LIBPR.0015 | Supersedes: Version 7 |
| Version: 8 | Page 10 of 17 |

Non Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

- 5.18 Make a 1:3 dilution of the reverse crosslinked sample.
- 5.19 Run the sample and the 1:3 dilution on an Agilent DNA 7500 assay to check fragment sizes. If no fragment is between 100-300bp, discuss with supervisor for further sonication.
- 5.20 Run a BR Qubit assay to determine the concentration.
- 5.21 Send the Agilent profile and Qubit quant to the supervisor to confirm there is enough chromatin to proceed with IPs. If there is not enough chromatin available to set-up IPs, proceed with Step 6.1-6.3 and then snap freeze the precleared chromatin in liquid nitrogen as per Appendix B.

6 Preclearing (Day 2)

- 6.1. Add 40µL of the blocked bead mix slurry to the pooled chromatin fraction, parafilm the tube and rotate at 4°C for 2 hrs.

Note: 40µL of blocked beads are sufficient for pre-clearing chromatin from 1x10e8 cells, adjust the amount of beads needed if different amount of cell number were used.
- 6.2. Spin at 4,000rpm for 2 min at 4°C. Transfer supernatant carefully into new 1.5mL non-stick tube labeled as precleared chromatin.
- 6.3. Spin the supernatant at 4,000rpm for 2 min at 4°C to ensure all beads have been removed.

7. Immunoprecipitation (IP) Set Up

- 7.1 Set up IPs in non-stick 0.6mL tubes. Use 5µg of chromatin/IP. Determine the volume of chromatin to use based on the Qubit concentration. Consult supervisor if not enough chromatin is present to set up the full panel of IPs or if the volume of 5µg of chromatin exceeds 30µL.
- 7.2 Fill out the ChIP Sample Summary Sheet found in :\\Library Core\Epigenomics\ChIP Sample Summary_121128-Present with the sample name, antibody details and Qubit quant of the chromatin to determine the amount of precleared chromatin and IP buffer + PIC to add to each IP.
- 7.3 Aliquot out the amount of precleared chromatin and IP buffer + PIC indicated on the ChIP Sample Summary Sheet.

Non Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

- 7.4 Add the corresponding antibody following the parameters listed on the table below (Table1). Parafilm the tubes and rotate at 4°C for 1 hr.

| Antibody Name | Company Name | Catalog # | ug or uL /IP |
|---------------|----------------|----------------------------|--------------|
| H3K4me1 | Diagnode | pAb-037-050 (C15410037) | 5ug |
| H3K4me3 | Cell Signaling | 9751S | 1ug |
| H3K9me3 | Diagenode | pAb-056-050 (C15410056) | 5ug |
| H3K27me3 | Diagenode | pAb-069-050 (C15410069) | 5ug |
| H3K36me3 | Abcam | ab9050 | 2.5ug |
| H3K27ac | Abcam | ab4729 | 2.5ug |

Table 1

Note: Antibody to be used in IPs might be different from those listed above. Refer to work request details from supervisor.

- 7.5 Pulse spin IP tubes.
- 7.6 Resuspend the Protein A/G prepared slurry by finger vortexing.
- 7.7 Add 20µL of Protein A/G prepared slurry to each IP using a wide bore P200 tip. Pipet up and down after adding beads to rinse the tip.
- 7.8 Parafilm tubes and rotate at 4°C overnight.
- 8. Washing (Day 3)**
- 8.1 Place ChIP Wash buffer and Final ChIP Wash buffer on ice.
- 8.2 Remove the parafilm and pulse spin the IPs.
- 8.3 Using a p200 tip with the end cut off, transfer the IPs to a Spin-X column.



| | |
|--------------------------------------|-----------------------|
| Chromatin Immunoprecipitation (ChIP) | |
| Document #: LIBPR.0015 | Supersedes: Version 7 |
| Version: 8 | Page 12 of 17 |

Non Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

- 8.4 Centrifuge the Spin-X columns at 1,000rcf for 1 min at 4°C. Discard the flow-through.
- 8.5 Add 400µL of ChIP Wash buffer to the Spin-X columns.
- 8.6 Rotate the beads at 4°C for 3 min.
- 8.7 Centrifuge the Spin-X column at 1,000rcf for 1 min at 4°C. Discard the flow-through.
- 8.8 Repeat Steps 8.5-8.7 with 400µL of ChIP Wash Buffer and then once with 400µL of Final ChIP Wash Buffer. Ensure there is no liquid on the beads after the last spin. If there is, repeat Step 8.7 until there is no liquid present on the beads.
- 8.9 Transfer the column to a new labeled 1.5mL non-stick tube.

9. Elution and Reverse Crosslinking

- 9.1 Pre-heat the thermomixer to 68°C.
- 9.2 Retrieve a 10uL aliquot of DTT from the -20°C freezer. Each aliquot is enough to make 2ml of the ChIP elution buffer + DTT mix. Scale up or down the number of DTT aliquots depending on the number of IPs set up. Store on ice until ready to add to the columns.
- 9.3 Add 2mL of ChIP Elution to each pre-aliquoted tube of DTT and warm up the solution on a thermomixer at 68°C.
- 9.4 Add 125 µL of pre-heated elution buffer + DTT to the columns.
- 9.5 Incubate the columns for 15 min at 68°C. Set the thermomixer to shake on intervals at 1,000rpm for 30 seconds, every 1 minute. Firmly attach the lid to the block. **DO NOT switch on the heat for the lid component.** Note the equipment ID of the thermomixer and the lid used.
- 9.6 Centrifuge the columns at 1,000 rcf for 1 min at room temperature. Ensure all liquid has spun through. Spin again if the liquid has not completely spun through the column. Do not change the collection tube.
- 9.7 Add another 125 µL of pre-heated ChIP elution buffer + DTT to the column. Incubate for 15 min at 68°C with shaking. Ensure the lid is on securely. **DO NOT switch on the heat for the lid component.**



| | |
|--------------------------------------|-----------------------|
| Chromatin Immunoprecipitation (ChIP) | |
| Document #: LIBPR.0015 | Supersedes: Version 7 |
| Version: 8 | Page 13 of 17 |

Non Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

- 9.8 Centrifuge the columns at 1,000 rcf for 1 min at room temperature. Ensure all liquid has spun through. Spin again if the liquid has not completely spun through the column. There should be 250 μ L of ChIP DNA. Detach the lid and discard the column.
- 9.9 Aliquot out 100ng for Input fraction into a non-stick 1.5mL tube. Store Input fraction at 4°C.
- 9.10 Add 250 μ L of elution buffer + DTT to the input sample (From Step 9.9).
- 9.11 Add 10 μ L of RNase A to all the IPs and input samples.
- 9.12 Invert tubes to mix and pulse spin.
- 9.13 Incubate the samples at 68°C on the thermomixer with the lid. **Switch 'ON' the heating component of the lid with NO SHAKING.** The lid is preset to 78°C. Ensure the lid is attached securely. Incubate the tubes for 5 hours or overnight.
- 9.14 Snap freeze the remaining precleared chromatin. See Appendix B.

10 Purification of DNA

- 10.1 Purify the DNA with Qiaquick columns using the QIAquick PCR Purification Kit. Use a 12.5mL Distriman pipettor for the distribution of the PBI and PE buffers.
- 10.2 Add 1250 μ L PBI to the 250 μ L sample and pipet up and down to mix. Add approximately 750 μ L max of the sample to the Qiaquick column.
- 10.3 Spin at 13,000rpm for 1 min. Discard the flow-through.
- 10.4 Add the remaining sample to the column.
- 10.5 Spin at 13,000rpm for 1 min.
- 10.6 Add 750 μ L of PE buffer.
- 10.7 Spin at 13,000rpm for 1 min. Discard the flow-through.
- 10.8 Spin column one additional time for 2 minutes at 13,000rpm.
- 10.9 Use p200 tip to aspirate any addition ethanol trapped on the rim in the Qiagen column.
- 10.10 Transfer the column to a new labeled 1.5mL non-stick tube.



| | |
|--------------------------------------|-----------------------|
| Chromatin Immunoprecipitation (ChIP) | |
| Document #: LIBPR.0015 | Supersedes: Version 7 |
| Version: 8 | Page 14 of 17 |

Non Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

- 10.11 Heat the EB buffer (Qiagen) to 50°C.
- 10.12 Let the column air dry for 1 minute with the lid open before adding EB.
- 10.13 Add 30 µL of EB to the centre of the column and let sit for 1 minute.
- 10.14 Spin down at 13,000rpm for 1 minute. Discard columns.
- 10.15 Wipe down work bench and small equipment with DNAway.
- 10.16 QC the ChIP samples by qPCR. Refer to LIBPR.0075 qPCR of ChIP Samples to set up qPCR.
- 10.17 Update the ChIP Sample Summary Sheet (see Step 7.2 for location).
- 10.18 Email the qPCR results to supervisors. Once the IPs and Input Control samples are approved for library construction proceed to Step 10.19.
- 10.19 Top up each IP and Input Control sample to 35µl with EB buffer.



| Chromatin Immunoprecipitation (ChIP) | |
|--------------------------------------|-----------------------|
| Document #: LIBPR.0015 | Supersedes: Version 7 |
| Version: 8 | Page 15 of 17 |

Non Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

Appendix A

| Reagent | Stock Concentration | Volume | Final Concentration |
|--------------------------|---------------------|-------------|---------------------|
| ChIP Lysis Buffer | | | |
| Tris-HCl (pH 8.0) | 1M | 25mL | 50mM |
| SDS | 20% | 25mL | 1% |
| EDTA (lab stock) | 0.5M | 10mL | 10mM |
| Ultrapure Water | | Up to 500mL | |
| Complete-Mini PIC tablet | | | Refer to protocol |

| | | | |
|--------------------------|------|-------------|-------------------|
| IP Buffer | | | |
| Tris-HCl (pH 8.0) | 1M | 5mL | 10mM |
| Triton X-100 | 100% | 5mL | 1% |
| Deoxycholate* | 10%* | 5mL | 0.1% |
| SDS | 20% | 2.5mL | 0.1% |
| NaCl | 5M | 9ml | 90mM |
| EDTA (lab stock) | 0.5M | 2mL | 2mM |
| Ultrapure Water | | Up to 500mL | |
| Complete-Mini PIC tablet | | | Refer to protocol |

| | | | |
|-------------------------|------|-------------|-------|
| ChIP Wash Buffer | | | |
| Tris-HCl (pH 8.0) | 1M | 10mL | 20mM |
| SDS | 20% | 2.5mL | 0.1% |
| Triton X-100 | 100% | 5mL | 1% |
| EDTA (lab stock) | 0.5M | 2mL | 2mM |
| NaCl | 5M | 15mL | 150mM |
| Ultrapure Water | | Up to 500mL | |

| | | | |
|-------------------------------|-------|-------------|-------|
| ChIP Final Wash Buffer | | | |
| Tris-HCl (pH 8.0) | 1M | 10mL | 20mM |
| SDS | 20% | 2.5mL | 0.1% |
| Triton X-100 | 100% | 5mL | 1% |
| EDTA (lab stock) | 0.5M | 2mL | 2mM |
| NaCl*** | 5M*** | 50mL | 500mM |
| Ultrapure Water | | Up to 500mL | |

| | | | |
|----------------------------|--------|-------------|-------|
| ChIP Elution Buffer | | | |
| NaHCO ₃ **** | 1M**** | 5mL | 100mM |
| SDS | 20% | 2.5mL | 1% |
| Ultrapure Water | | Up to 500mL | |

*10% Deoxycholate stock Make a 10% stock by adding 1g of powder and top up to 10mL with water in a 50mL falcon tube



| | |
|--------------------------------------|-----------------------|
| Chromatin Immunoprecipitation (ChIP) | |
| Document #: LIBPR.0015 | Supersedes: Version 7 |
| Version: 8 | Page 16 of 17 |

Non Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

***NaHCO₃ Make a 1M stock by adding 8.4g of powder and top up to 100mL with water in sterile bottle

OFFICIAL



| | |
|--------------------------------------|-----------------------|
| Chromatin Immunoprecipitation (ChIP) | |
| Document #: LIBPR.0015 | Supersedes: Version 7 |
| Version: 8 | Page 17 of 17 |

Non Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

Appendix B

Snap Freezing Precleared Chromatin

1. Record the concentration of the precleared chromatin from the qubit quant in Step 5.18 on to the precleared chromatin tube.
2. Measure the remaining volume of precleared chromatin and enter the volume in the ChIP Sample Summary Sheet (File location in Step 7.2).
3. Put the tubes into an empty pipette tip box and transport the tubes to the 5th floor tissue culture room.
4. With proper personal protective gear, including face shield, remove one of the cell storage chambers from the liquid nitrogen cryo tank.
5. Place the tubes into the chamber and lower the chamber into the liquid nitrogen cryo tank.
6. The tubes should snap freeze immediately. If unsure, leave the tubes in the tank for a minimum of 1 minute.
7. Transfer the tubes to an appropriate box in the -80°C located on the 5th floor. Record the rack location and enter it into the ChIP Sample Summary Sheet.



| | |
|--------------------------------------|-----------------------|
| Chromatin Immunoprecipitation (ChIP) | |
| Document #: LIBPR.0015 | Supersedes: Version 7 |
| Version: 8 | Page 18 of 17 |

Non Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

Appendix C

LIMS for ChIP

1. SLX-ChIP: CHP_version # pipeline (eg. CHP_2.1)
2. Update antibody volume for control in LIMS under the “Antibody_Inventory_Tracking” view.

OFFICIAL