

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 1 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

Native ChIP Using 100,000 Cells

I. Purpose

To provide details for immunoprecipitation reactions using 100,000 cells from cells or tissues.

II. Scope

All procedures are applicable to the BCGSC Library Core and Library TechD groups.

III. Policy

This procedure will be controlled under the policies of the Genome Sciences Centre, as outlined in the Genome Sciences Centre High Throughput Production Quality Manual (QM.0001). Do not copy or alter this document. To obtain a copy see a QS associate.

IV. Responsibility

It is the responsibility of all personnel performing this procedure to follow the current protocol. It is the responsibility of the Production Coordinator to ensure personnel are trained in all aspects of this protocol. It is the responsibility of Quality Systems to audit this procedure for compliance and maintain control of this procedure.

V. References

Document Title	Document Number
ndChIP-seq	EGL017.02
Preparation and QC of SeraMag Bead Solution	EGL008.2

VI. Related Documents

Document Title	Document Number
N/A	

VII. Safety

All Laboratory Safety procedures will be complied with during this procedure. The required personal protective equipment includes a laboratory coat and gloves. See the material safety data sheet (MSDS) for additional information.

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 2 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

VIII. Materials and Equipment

Name	Supplier	Number	Model or Catalogue #	
Small size safetouch nitrile gloves	Ultident	296359683		✓
wet ice	In house	N/A	N/A	N/A
Ice bucket	Fisher	11-675-58		✓
1.5 mL Microtubes	Diamed	PRE150-B		✓
15 mL Conical Tubes	BD Falcon	352097		✓
50 mL Conical Tubes	BD Falcon	352070		
Gilson P10 pipetman	Mandel	GF-44802		✓
Gilson P20 pipetman	Mandel	GF23600		✓
Gilson P200 pipetman	Mandel	GF-23601		✓
Gilson P1000 pipetman	Mandel	GF-23602		✓
Diamond Filter tips DFL10 (10 Tipacks of 96 racked filter tips)	Mandel Scientific	GF-F171203		✓
Diamond Filter tips DFL30 (10 Tipacks of 96 racked filter tips)	Mandel Scientific	GF-F171303		✓
Diamond Filter tips DFL200 (10 Tipacks of 96 racked filter tips)	Mandel Scientific	GF-F171503		✓
Diamond Filter tips DFL1000 (10 Tipacks of 96 racked filter tips)	Mandel Scientific	GF-F171703		✓
Galaxy mini-centrifuge	VWR	37000-700		✓
Large Kimwipes (Kimberly Clark/Kimtech)	Fisher Scientific	06-666-1A		✓
Black ink permanent marker pen	VWR	52877-310		✓
Small Autoclave waste bags 10"X15"	Fisher Scientific	01-826-4		✓
Ultra Pure Water (RNase/DNase free)	Invitrogen	10977-023		✓
Anhydrous Ethyl Alcohol (100% Ethanol)	Commercial Alcohols	People Soft ID: 23878		✓
DNA away	Molecular Bioproducts	7010		✓
1M Tris Solution, pH 8.0	Ambion	AM9856		✓
Triton® X-100, laboratory grade	Sigma	X100-100ML		✓
20% SDS Solution	Ambion	AM9820		✓
Centrifuge, Eppendorf 5417R, refrigerated high-speed, 115V	Fisher Scientific	5417 R	✓	
VX-100 Vortex Mixer	Rose Scientific	S-0100	✓	
Eppendorf Thermomixer 1.5 mL	Eppendorf	21516-166	✓	
Parafilm	Fisher Scientific	13-374-12		✓
Water Bath	Fisher Scientific	Isotempn 220	✓	
5M NaCl	Ambion	AM9760G		✓
1M DTT	Invitrogen	P2325		✓
1M Tris HCl pH 7.5	Invitrogen	15567-027		✓
0.5M EDTA	Ambion	AM9260G		✓
Deoxycholic Acid, Sodium Salt	Fisher Scientific	AC218590250		✓
Sodium Bicarbonate	Sigma Aldrich	S5761-500G		✓
PCR Clean ChIPseq Beads with 30% PEG	Aline Biosciences	C-1007		✓

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

0.2 mL Ultra Rigid Skirted 96 well PCR plate	Thermoscientific	FSSP9741450		✓
Micrococcal Nuclease	NEB	M0247S		✓
Nonstick 1.5 mL tubes	Ambion	AM12450		✓
1M DTT	Invitrogen	P2325		✓
RNase Away	Invitrogen	10328-011		✓
Dynabeads Protein G, 5 mL size	Invitrogen	100-04D		✓
Dynabeads Protein A, 5 mL size	Invitrogen	100-02D		✓
Protease Inhibitor Cocktail	Calbiochem	539134		✓
LabQuake Shaker/Rotator with clips	Barnstead	415110	✓	
Buffer G2	Qiagen	1014636		✓
Qiagen Protease	Qiagen	19155		✓
DynaMag 2 Magnet	Invitrogen	12321D		✓
Disposable Trough	VWR	21007-972		✓
Multi 12-channel Pipette P20	Rainin	17013803		✓
Multi 12-channel Pipette P200	Rainin	17013805		✓
EB Buffer	Qiagen	19086		✓
96-well Plate Magnet	Alpaqua	002523		✓
Tape Pads: Adhesive Plate Sealer	Qiagen	19570		✓
Aluminum Foils	VWR	60941-126		✓
Rainin 200µl barrier tips	Rainin	RT-L200F		✓
Rainin 20µl barrier Tips	Rainin	RT-L20F		✓
Sodium Butyrate	Millipore	19-137		✓
Sodium Azide	Sigma-Aldrich	S2002-100G		✓
Liquid Nitrogen	Praxair	NI M-FILTER		✓
Dry Ice	In house	N/A		✓
Steel plate	In house	N/A		✓
Liquid Nitrogen Cooled Mini Mortar and Pestle Set, Bel-Art	VWR	89233-994		✓
Single Edge Industrial Razor Blades	VWR	55411-050		✓
Sterile Disposable Forceps	VWR	12576-934		✓
Sterile Cell Culture Dishes (150mmx25mm)	Fisher Scientific	08-772-25		✓
Sterile Cell Culture Dishes (60mmx16mm)	Fisher Scientific	08-772-31		✓
H3K4me1 Rb Antibody	Consult Supervisor	N/A		
H3K4me3(C42D8) Rb Antibody	Consult Supervisor	N/A		
H3K9me3 Rb Antibody	Consult Supervisor	N/A		
H3K27me3 Rb Antibody	Consult Supervisor	N/A		
H3K36me3 Rb Antibody	Consult Supervisor	N/A		
H3K27ac Rb Antibody	Consult Supervisor	N/A		
PCR tube strip, Domed 12 cap strip	Biorad	TCS 1201		✓
Eppendorf Thermomixer C	Fisher Scientific	05-412-503		✓
96 Well PCR plate Smartblock for Thermomixer C	Fisher Scientific	05-412-512		✓
T36 Disinfect	VWR	CA11007-034		✓
PEG-8000 (molecular biology grade)	Sigma-Aldrich	P5413-500G		✓
Seramag speed beads	Fisher Scientific	09-981-123		✓
Portable Pipet Aid, Multispeed XP, rechargeable	Fisher Scientific	13-681-15E		✓

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 4 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

50 mL conical tube	VWR	CA21008-940		✓
15 mL conical tube	VWR	CA21008-918		✓

IX. Procedure

1. Retrieval of reagents and equipment preparation

- 1.1. Put on a lab coat and clean pair of gloves.
- 1.2. Wipe down the assigned specific workstation, pipettes, and small equipment.
- 1.3. Lay down new benchcoat.
- 1.4. Change gloves.
- 1.5. Confirm with supervisor as to sample type (tissue or cells).
 - 1.5.1. If sample source is cells, verify the number of cells/vial.
- 1.6. Prepare worksheets keeping in mind:
 - What is the sample source (i.e. tissue or cells)
 - Number of IPs and what are the targets/antibodies to IP
 - Include positive control IPs (100K HL60 cells:whole histone modified panel of antibodies)
- 1.7. Note that this protocol is optimized for **100,000 cells (100K)/IP**.
- 1.8. A minimum of 2L of liquid nitrogen is required for tissue samples.
- 1.9. A small amount of dry ice is required for tissue samples.
- 1.10. Dilute antibodies if needed. See Appendix B if dilutions are required. Can be prepared ahead of time and stored at -20°C.
- 1.11. Prepare ahead of time 30% PEG SeraMag beads as detailed in Appendix C.

Ensure all buffers, reagents and enzymes have been validated using the Native ChIP production pipeline up to library qPCR QC.
- 1.12. Prepare buffers and reagents and keep on ice

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 5 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

DAY 1

2. Preparation of Ab-Bead Complex

Prior to starting protocol ensure the following:

Prepare worksheets and plate layout before starting protocol

Ensure there is enough diluted stock of antibodies. See Appendix B if not

If probing with H3K27ac, add 10mM (final) of sodium butyrate to separate aliquots of IP buffer + PIC and Lysis Buffer + PIC. Use this buffer specifically for H3K27ac IP wells as sodium butyrate is a deacetylation inhibitor.

- 2.1. Add 20µL of Calbiochem protease inhibitor cocktail (500x) into 10mL of IP Buffer. (diluting PIC to 1x). Mix well and keep on ice.
- 2.2. Retrieve Dynabeads Protein A and Dynabeads Protein G from the 4°C and mix each bottle very well to achieve a homogenous bead slurry. Can be gently vortexed. A 50:50 bead slurry (volume of Dynabeads Protein A:volume of Dynabeads Protein G) is required.
- 2.3. Label a 1.5mL non-stick tube with '50:50 A/G beads' and the date.
- 2.4. Prepare a 50:50 bead slurry of Protein A:Protein G by adding equal amount of Dynabeads Protein A and G.
- 2.5. Transfer appropriate amount of beads to the labeled 1.5mL non-stick tube by using the following formula:

$$(23\mu\text{L Protein A} + 23\mu\text{L Protein G}) \times \# \text{ of IPs} + \text{dead volume}^* = \text{amount of beads}$$

*dead volume is 3 reactions worth of beads
- 2.6. Place bead tube on magnet stand and let separate. Remove supernatant.
- 2.7. Remove bead tube from magnet stand and place on ice.

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

- 2.8. Add equal volume of IP Buffer + PIC to beads. (1:1 v/v bead:IP buffer)
- 2.9. Mix by pipetting up and down. Ensure thoroughly mixed.
- 2.10. Place bead tube back on magnet stand and let separate. Remove supernatant.
- 2.11. Repeat wash 2 more times with IP + PIC buffer.
- 2.12. Resuspend beads in equal volume of IP Buffer + PIC mix. Keep on ice.
- 2.13. On ice, aliquot 46µL of washed Dynabeads per well of sample into a new plate.
- 2.14. On another new plate, aliquot 130µL per well of IP buffer + PIC. Can use a reagent trough for this.
- 2.15. Add 20µL of washed Dynabeads to each well containing the 130µL of IP buffer + PIC. Mix by pipetting up and down.
- 2.16. Seal remaining Dynabeads in plate and keep on ice. These will be used for pre-clearing chromatin later.
- 2.17. Addition of Antibodies:

If necessary dilute the stock antibody to a concentration that is reasonable to pipette. See Appendix B if dilution of stock antibody is necessary. Do not pipette less than 1µL. Add 0.5-0.75µg of antibody to the wells containing Dynabeads and IP buffer + PIC (plate from Step 2.15) to generate the Ab-bead complex. See Table 1 below for the amount of antibody to add. Ensure that a positive control IP has been set up as well using H3K4me3 antibody.

Name	Amount (µg)
H3K4me3	0.75
H3K4me1	0.5
H3K27me3	0.5
H3K36me3	0.5
H3K27ac	0.5
H3K9me3	0.5

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 7 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

Table 1: Amount of Antibody

- 2.18. Set multi-channel P200 pipette to 100µL and mix each row up and down 10x. Change tips between each row.
- 2.19. Cover the Ab-bead complex plate with BioRad domed cap strip lids securely. Ensure the beads are mixing thoroughly and no air pockets are visible in the wells. If there are air pockets, gently flick the lid(s) of the wells(s) to disrupt the air pocket.
- 2.20. Incubate the plate at 4°C on a rotating platform for a minimum of 2.5 hrs or until pre-clearing is done.
- 2.21. Proceed to Step 5 after incubation.

3. Cell Lysis and MNase I Digestion of Chromatin

- 3.1. Retrieve lysis buffer from 4°C and place on ice.
- 3.2. Aliquot 10mL into 15mL tube.
- 3.3. Add 20µL of Calbiochem Protease Inhibitor cocktail (500x) to the buffer. Mix and keep on ice.
- 3.4. If working with H3K27ac antibody, aliquot 990µL of lysis buffer + PIC into 1.5mL tube. Add 10µL of 1M sodium butyrate (final 10mM). Mix well and put on ice.
- 3.5. Prepare sample for lysis and MNase I digestion. **For fresh or frozen 100K cells follow Steps 3.6. For tissue samples follow Step 3.7.**

3.6. For Fresh/Frozen Cells:

- 3.6.1. Verify the number of cells/vial prior to proceeding as mentioned in Step 1.5.1. Consult with supervisor if unsure.
- 3.6.2. Depending on the number of cells/vial and the number of IPs to set up:

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 8 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

Thaw 100K cells /1 IP + 1 extra IP (to account for loss) at 37°C for a few seconds in a water bath.

****For a complete histone panel, use 700,000 cells for 100K IPs. Equally divide the extra 100K cells amongst the 6 IPs of a whole panel; 117K cells / well****

*****Any reference to a '100K' IP using a cells sample includes the extra 17K cells*****

Note: If the sample has more than the # of cells required for the Native ChIP process, supervisor will advise to split the frozen cell stock. In that case, please refer to the following instruction in step 3.6.2.1-3.6.2.13:

- 3.6.2.1. In a 15 ml falcon tube, pipette 5mL of PBS and add 10uL of 500x PIC. Keep on ice.
- 3.6.2.2. Prechill Eppendorf centrifuge to 4°C.
- 3.6.2.3. Retrieve the source cell tube in a small styrofoam box containing dry ice.
- 3.6.2.4. Label and prechill 1.5mL eppendorf tubes for each cell aliquot.
- 3.6.2.5. Thaw the source tube quickly either by hand or by dipping in a 37°C water bath.
- 3.6.2.6. Gently flick the source tube to ensure all the cells are thawed and evenly suspended.
- 3.6.2.7. Add an appropriate amount of PBS containing 1xPIC based on the number of cells reported by the collaborator.
- 3.6.2.8. Mix the cells by pipetting up and down gently to ensure there are no cell clumps.
- 3.6.2.9. Aliquot out, the volume of cells needed into a previously prechilled 1.5ml eppendorf tube.

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 9 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

- 3.6.2.10. Spin down both the source and cell aliquot tubes at 2655g (5,000rpm in Eppendorf Centrifuge 5417R), for 5min at 4°C.
- 3.6.2.11. Pipette off and discard the supernatant from both the source and aliquot cell tubes.
- 3.6.2.12. Place the source tube back on dry ice.
- 3.6.2.13. To the cell aliquot, proceed to Step 3.6.3.
- 3.6.3. Top up each 100K worth of cells up to 30µL with cold Lysis buffer + PIC. Scale up if necessary. Pipette up and down at least 10 times to mix. Make sure there are no clumps. Remember to add lysis + PIC+10mM sodium butyrate to pellet for H3K27ac IP.
- 3.6.4. Also thaw out a pellet of HL60 cells as a positive control. Follow Step 3.7.20. Add cells accordingly as set up on plate layout worksheet.
- 3.6.5. Transfer each 100K cell lysate to 1 well each of a 96 well plate. Refer to the prepared sample table (See worksheet) to ensure the correct lysate is transferred to the Ab-bead complex plate.
- 3.6.6. Cover the plate with a plastic seal and incubate for 20 min. During incubation prepare MNase I dilution and master mix, Step 3.8.
- 3.7. For Tissue Samples:**
- 3.7.1. Ensure 2L of liquid nitrogen and a small quantity of dry ice is present in the lab before starting.
- 3.7.2. To cut a piece of tissue from the original source, collect the following items per 1 **sample** and put them in the -80°C for at least 1 hour prior to cutting the tissue. It is also acceptable to put the apparatus and supplies in the -80°C overnight. **It is**

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 10 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

crucial that all reagents and supplies used in cutting and grinding of the tissue remain extremely cold.

- Steel plate
- 2 razor blades
- 1 disposable forcep
- 1 small petri dish
- large lid of a styrofoam box
- 1 non-stick 1.5mL tube
- 15mL tube (to hold disinfected pestle)
- Mini Mortar and Metal Pestle Set (ensure it is thoroughly disinfected with DNase away, rinsed with sterile ultrapure water, and then wiped with 70% ethanol.

***If working with multiple tissue samples, gather all the supplies needed. See Steps 3.7.23-3.7.30 for additional set up of processing multiple tissue samples.**

- 3.7.3. Prepare the BSC by turning it on for 15-20min. After purging the BSC, wipe it clean with T36 disinfectant. Wipe with 70% ethanol to remove traces of T36 disinfectant. Tape a small biohazard waste bag to the inside wall to dispose of biohazard waste.
- 3.7.4. Wipe down the bottom of the non-consumable supplies with 70% ethanol and transfer all items from Step 3.7.2 into the BSC carefully. It is not necessary to wipe down the steel plate as it was previously decontaminated.
- 3.7.5. Transfer dry ice into the large styrofoam lid, in a single layer. Lay the steel plate onto the dish. Onto the lid containing dry ice, place the remaining supplies except the mini mortar. Ensure the supplies are fully covered by dry ice.
- 3.7.6. Carefully transfer liquid nitrogen into the mortar bowl. Fill the bowl about half way to avoid spilling once the mortar is placed back onto it. Tightly secure the mortar onto the bowl to prevent rapid evaporation of the liquid nitrogen. Place the

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 11 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

1.5mL non-stick tube in the tube holder of the mortar bowl. See Figure 1 for set up.



Figure 1

- 3.7.7. Wait until the mortar and the bowl become 'frosty' to start. This will take around 5 minutes. The apparatus will stay extremely cold for about 10 minutes before the liquid nitrogen needs to be topped up. It is important to work fast but safely.
- 3.7.8. Retrieve the source tissue on dry ice. Leave on dry ice until ready to cut.
- 3.7.9. Place the small petri dish onto the steel plate, and leave open. Carefully remove paper protection wrap from 1 or 2 razor blades and place in the lid of the petri dish. Try to minimize handling to prevent warming up of supplies.
- 3.7.10. Change gloves and double gloving is optional. Transfer the source tissue onto the small petri dish. Cut a small piece using the chilled razor blade(s), approximately 2mm^3 . Be extremely careful when cutting the tissue. If necessary, use the lid of the petri dish as a shield to deflect any tissue(s) piece(s) that become airborne. 2mm^3 of tissue is required only if setting up a whole histone modification panel (6 IPs). Adjust the amount of tissue to be cut if setting up fewer or more IPs. Quickly transfer the source tissue back into source tube and put on dry ice. Consult with supervisor if unsure of size amount.
- 3.7.11. Transfer the cut tissue into the chilled 1.5mL non-stick tube sitting in the mortar bowl using the razor blade or disposable forceps. Leave the lid open.
- 3.7.12. Using the chilled metal pestle, apply controlled yet substantial pressure against the frozen tissue piece while turning the pestle within tube. The tissue should pulverize into a powder within a few turns.

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 12 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

- 3.7.13. **DO NOT** remove the tube from the tube holder of the mortar as this action will thaw the tissue. If it is absolutely necessary to check the grinding efficiency, place the tube immediately on dry ice. Also keeping the metal plunger in the air will warm it up as well. Be aware of all supplies and apparatus as this process is only successful if everything is very cold, at the same time ensure your safety.
- 3.7.14. After the grinding, place the tube on dry ice.
- 3.7.15. Remove the sample tube from dry ice and carefully add 30 μ L of Lysis Buffer + PIC/IP. Scale up if required as noted below. Once the ground up tissue is resuspended, it is safe to put the tube on wet ice. With the tube in one hand, ensure the ground up tissue thaws as it is being lysed at the same time.
- *For an entire histone modification panel: (6 IPs + 1 extra for loss)
30 μ L=210 μ L Lysis Buffer + PIC***
- 3.7.16. Pulse spin very gently if necessary to bring down any lysis buffer from the sides of the tube. Leave sample tube on wet ice for 10min once the ground up tissue has thawed. During the 10min incubation proceed to Step 3.7.20 at the 6 min mark.
- 3.7.17. Pulse spin the tube on a microfuge to pull down cellular debris. Gently put back on ice to avoid disturbing the pellet.
- 3.7.18. Carefully and thoroughly transfer the supernatant, which contains the chromatin, to a fresh 1.5mL non-stick tube. Avoid disturbing the pellet.
- 3.7.19. Put sample tube back on wet ice and equally distribute the supernatant amongst wells assigned for an IP of a 96-well plate. See plate layout worksheet. For example, if the entire histone modification panel of antibodies is going to be used then equally divide the supernatant amongst 6 assigned wells.
- 3.7.20. Pull out HL60 control cells from the -80°C freezer on dry ice. Each cell pellet is at 1 million cells/100 μ L of lysis buffer +PIC. Add 150 μ L of lysis buffer + PIC, thaw and resuspend cells simultaneously. Submerge the tube in a 37°C water bath for a quick thaw, is an option as well. Aliquot 30 μ L of cell suspension to appropriate assigned wells, according the plate layout worksheet. There will be approximately 117K cells/well for a 100K IP.

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 13 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

3.7.21. Cover the plate with a seal and let incubate on ice for 20min. It is safe to go over by 5-10 minutes since there is PIC present in the lysis buffer.

3.7.22. During the 20min incubation, proceed to Step 3.8-3.16 BUT DO NOT ADD the micrococcal nuclease until ready to proceed to the digestion step. After the 20min incubation is complete, proceed to **Step 3.17**.

3.7.23. If processing multiple tissue samples it is necessary to keep the pulverized tissue frozen as subsequent samples are being cut and crushed. Until all samples are processed can they be thawed all at once for lysis and subsequent steps. Along with the supplies listed in Step 3.7.2, a few additional items are required and are listed below:

- small container of liquid nitrogen close by in case a top up is required
- 15mL falcon tube(s) containing 9mL DNase away
- 1 falcon tube to hold disinfected pestle
- squirt bottle with ultrapure water
- P200 pipette and filtered tips
- lid of a large styrofoam box filled with 1 layer of dry ice pellets
- small styrofoam box filled with wet ice to contain lysis buffer
- 1 non-stick tube/sample pre-labeled with sample ID on lid and side of tube with a permanent marker

3.7.24. See Figure 2 below as an example of processing multiple tissue samples. As visualized a substantial amount of supplies are required. Ensure the set up follows proper aseptic techniques as best as possible and the flow of work starts at one end of the BSC and ends at the other, as directed in Figure 2 as well. This also provides guidelines to ensure safety of the user as well. This allows minimum disruption of the airflow as possible.



Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 14 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*




Figure 2: Additional Set Up For Processing Multiple Tissue Samples

- 3.7.25. Retrieve the tissue source tubes on dry ice. Process 1 sample at a time. After cutting the tissue sample, place any remaining piece(s) back into source tube. Immediately put back on dry ice to avoid thawing. Once the cut tissue piece has been pulverized, snap the lid shut while the tube is still in the mortar.
- 3.7.26. By holding the tube firmly on the top (closest to the lid), transfer it to box containing dry ice. This is to prevent the tissue from warming up.
- 3.7.27. Immerse the metal portion of the pestle used to grind the tissue piece into a tube containing DNase away. Leave for about 1 minute to disinfect. Pull it out, rinse the pestle with ultrapure water and wipe dry with a kimwipe. Remove traces of water by spraying with 70% ethanol and wipe dry. Put the pestle in a clean 15mL tube that has been sitting in dry ice to keep thoroughly cold.
- 3.7.28. Change gloves and proceed to the next sample.
- 3.7.29. Once all samples have been processed, safely remove all tubes from liquid nitrogen and add 30µL of lysis buffer + PIC/IP. Scale up as required. Resuspend slowly by pipetting up and down and place on wet ice to let completely thaw. Continue protocol with Step 3.7.16.
- 3.7.30. When ready, carefully and thoroughly clean and tidy the BSC, discarding all disposable supplies that came into contact with tissue into the biohazard waste

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 15 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

bag. Dispose of razors into the sharps container. Soak the mortar bowl, pestle and steel plate in DNase Away and then in ultra pure water to disinfect. Remove traces of water by spraying with 70% ethanol.

- 3.8. Prepare the MNase I dilution buffer in a 1.5mL tube, as according to Table 2. Make sure to mix it very well and keep it on ice. Make the buffer **fresh** for every use.

Reagent	Amount (μL)
1M Tris, pH 7.5	10
5M NaCl	10
0.5M EDTA	2
Glycerol	500
Ultrapure Water	478

Table 2: MNase I dilution buffer

- 3.9. Dilute Stock MNase enzyme (2000U/μL):

Amount of Cells	Dilution Factor	Amt Dilution Buffer (μL)	Amt Stock Enzyme (μL)
100K	1/10 (200U/μL)	18	2

Table 3: Dilution of Stock MNase Enzyme

- 3.10. Mix well by pipetting gently up and down. Scale up if needed. Avoid pipetting less than 2μL of enzyme. Keep on ice.

- 3.11. Prepare the MNase I digestion Master Mix in a 1.5mL tube, based on the worksheet.

Reagent	Volume (μL)/Rxn
20mM DTT*	1

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 16 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

MNase (200U/μL)	0.45
10X MNase Buffer	4
Ultrapure Water	4.5
Total/Reaction	10

Table 4: MNase I Digestion Master Mix

*1M DTT stock is diluted to 200mM and then to 20mM with water

- 3.12. Mix the brew very well by pipetting up and down slowly. Keep on ice.
- 3.13. Turn on thermoblock and set temperature to 25°C.
- 3.14. In a new 96 plate aliquot 10μL of the MNase I Master Mix per each row of samples plus 5μL of dead volume. Keep on ice.
- Example: For 2 rows the volume should be $(10\mu\text{L} \times 2) + 5\mu\text{L} = 25\mu\text{L}/\text{well} \times \# \text{ of IPs}$.
- 3.15. Aliquot 1% triton/1% DOC into a new 96-well plate. Label the plate. Aliquot 5.5μL per each row of samples plus 5μL of dead volume. Keep on ice.
- Example: For 2 rows the volume should be $(5.5\mu\text{L} \times 2) + 5\mu\text{L} = 16\mu\text{L}/\text{well} \times \# \text{ of IPs}$.
- 3.16. In a new 96 well plate, Aliquot 5.5μL of 250μM EDTA per each row of samples plus 5μL of dead volume. Keep at room temperature.
- Example: For 2 rows the volume should be $(5.5\mu\text{L} \times 2) + 5\mu\text{L} = 16\mu\text{L}/\text{well} \times \# \text{ of IPs}$.
- 3.17. After the lysates finish incubating for 20min on ice, remove the plate from ice and put on bench.
- 3.18. It is important to ensure the plate is at room temperature before proceeding to the next step.
- 3.19. If necessary, pulse spin the plate at room temp at 1500rpm.
- 3.20. Place the plate on the thermoblock preheated to 25°C. Using a P20 multichannel pipette, add 10μL of MNase I Master Mix to the first row of samples.

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 17 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

3.21. Mix by pipetting up and down 15x. Start the timer for 6 min. Gently place a cover seal on the plate and let incubate at 25°C on a thermoblock with lid on.

3.22. Proceed to the second row, if required. Make note of the mixing time.

Note: Do not process more than 2 rows at a time to make sure that each digestion is exactly for 6 min.

3.23. At the halfway mark of the incubation time, remove seal and mix by pipetting up and down 10x. Cover with a plastic seal. Change tips between rows.

3.24. Add 5.5µL of 250µM EDTA, with the plate on the thermoblock, to stop the digestion using a P20 multichannel pipette. Set the pipette to 20µL and mix 15x by pipetting up and down at room temp. Repeat with the second row if necessary. Ensure timing is equal between rows.

Note: This step should be performed so that each row of samples is digested with MNase I for exactly 6 min.

3.25. Repeat the process of MNase I digestion + EDTA quenching 2 rows at a time if processing more than 2 rows of samples.

3.26. When completed place the plate on ice.

3.27. Using a P20 multichannel, add 5.5µL of 1% Triton/1% DOC into each row of samples. Adjust the pipette volume to setting to 20µL and mix 15x by pipetting up and down.

3.28. Seal the plate with a plastic cover, and spin at 200g for 1min at 4°C.

3.29. Incubate on ice for 20min. Proceed to Step 4.1 during incubation time.

4. Input Separation and Preclearing

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 18 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

- 4.1. While the MNase I digestion plate is incubating, aliquot 130µL of IP Buffer + PIC into 1 row of a new 96-well plate, if working with just one row. Seal plate with plastic cover and keep on ice. If processing multiple rows of samples pour the IP Buffer + PIC into a disposable plastic reservoir. (Note: 130µL IP Buffer + PIC + Sodium Butyrate for H3K27ac IP).
- 4.2. After 20min incubation, using a P200 pipette pool all samples (of the same cell type) into a non-stick 1.5mL tube. Mix well by gently inverting tube and pulse spin. Transfer 12µL of the pool into a new non-stick 1.5µL tube. This will be the INPUT DNA for the sample. Store at 4°C until the following day.
- 4.3. Using a P200 pipette measure the volume of each pool. Note the volume. Distribute the sample equally (amongst number of IPs to set up) into a new 96 well plate. Put plate on ice. (See worksheet for plate layout)
- 4.4. Add 100-110µL of IP Buffer + PIC (from Step 4.1) to each sample well and mix 10x by pipetting up and down. If the volume of the chromatin/IP is high (>50 µL) add 100 µL. Add IP Buffer+PIC+Sodium Butyrate for H3K27ac IPs. Make sure that each well is well mixed before proceeding to the next step. Keep plate on ice.
- 4.5. Retrieve the pre-washed and aliquoted Dynabeads that were kept on ice from Step 2.16.
- 4.6. Using a P20 multichannel pipette, mix the beads by pipetting up and down a few times. Try not to introduce bubbles.
- 4.7. Using a P20 multichannel pipette, add 20µL of beads to each well of sample. After each addition, mix up and down a few times. Use fresh tips for each row.
- 4.8. To make sure that the samples are well mixed, use a P200 multichannel pipette set to 100µL and mix each row 10x by pipetting up and down. Change tips between rows.
- 4.9. Seal the plate very well with Biorad domed cap strip. Ensure there are no air bubbles present in the wells. If there are, gently flick the well to pop it.

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 19 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

4.10. Incubate on a rotating platform at 4°C for 1.5hr.

5. Immunoprecipitation Reaction

- 5.1. Place the Ab-Bead complex (from Step 2.20) on a magnet and wait for the beads to separate.
- 5.2. Carefully remove and discard the supernatant.
- 5.3. Remove the plate from the magnet and put the plate on ice.
- 5.4. Place the pre-clearing reaction plate (from Step 4.10) on a magnet and wait for the beads to separate.
- 5.5. Set the P200 multichannel pipette to 160µL and transfer the supernatant from the pre-clearing reaction plate to the Ab-bead complex plate kept on ice. Ensure that no beads are transferred.
- 5.6. Repeat the transfer and use a low volume multi channel pipette to ensure all of the supernatant from the pre-clearing reaction plate has been transferred to the Ab-bead complex plate.
- 5.7. Gently mix 15x by pipetting up and down. Ensure no air bubbles are present in the wells. Gently flick the plate if there are to pop them.
- 5.8. Seal the plate well with domed cap strip and incubate at 4°C on a rotating platform over night.

DAY 2

6. Washes and Elution

- 6.1. Put on a clean lab coat, fresh gloves, wipe down bench and pipettes with RNase Away solution. Turn on and set thermomixer temperature setting to 65°C.
- 6.2. Change gloves.

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 20 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

- 6.3. Retrieve ice and place Low Salt Wash Buffer and High Salt Buffer on ice.
- 6.4. Prepare Elution Buffer as described in Appendix A.
- 6.5. Aliquot 40µL of Elution Buffer per row of samples into a new 96-well plate. Cover and store at RT.
- 6.6. Spin down the IP reaction plate (from Step 5.8) at 200g for 1min at 4°C.
- 6.7. Place the IP reaction plate on a magnet and wait for the beads to separate.
- 6.8. Carefully remove and discard the supernatant.
- 6.9. Remove the IP reaction plate from the magnet and place on ice.
- 6.10. Pour out 3mL of Low Salt Wash buffer per row of samples into a trough and keep on ice.
- 6.11. Add 120µL of Low Salt Wash buffer per well to the IP reaction plate.
- 6.12. Slowly mix 10x to resuspend the beads.
- 6.13. Place the IP reaction plate on a magnet.
- 6.14. Set the multichannel pipette at 130µL and discard the supernatant.
- 6.15. Repeat the Low Salt buffer wash for a total of 2 washes.
- 6.16. Place the IP reaction plate on ice.
- 6.17. Pour out 3mL of High Salt Wash buffer per row of samples into a trough and keep on ice.
- 6.18. Add 120µL of High Salt Wash buffer per well to the IP reaction plate.
- 6.19. Slowly mix 10x to resuspend the beads.

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 21 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

- 6.20. Place the IP reaction plate on a magnet.
- 6.21. Set the multichannel pipette at 130µL and discard the supernatant.
- 6.22. Remove the IP reaction plate from the magnet and place it on ice. Place a new plate next to it.
- 6.23. Add 120µL of High Salt Wash buffer to each well of the IP reaction plate.
- 6.24. Slowly mix 10x to resuspend the beads.
- 6.25. Transfer the resuspended beads into the new plate. Place on magnet. Do not discard the source plate yet, keep on ice.
- 6.26. Once the bead suspension has started to clear, aspirate 20µL of supernatant from the new plate in Step 6.25 and rinse out the wells of the source plate to make sure all the beads are transferred. Ensure the aspirant is going into the corresponding wells from new plate to source plate and vice versa. Discard source plate.
- 6.27. Aspirate supernatant from new plate discard the supernatant.
- 6.28. Take the plate off the magnet and keep the plate at room temperature.
- 6.29. Add 30µL of Elution buffer to each well and mix slowly 10x to resuspend the beads.
- 6.30. Seal the plate with domed cap strip and incubate at 65°C for 1.5hrs with a mixing speed of 1350rpm.
- 6.31. Check the elution reactions every 15 minutes for the first 30min to ensure the beads do not settle. In this step we want to agitate the beads enough to prevent them from settling while avoiding splashing up of the liquid. If the beads are on the bottom, increase the mixing step to 1800rpm for a few seconds.
- 6.32. Spin down the plate at 200g for 1 min at RT.
- 6.33. Place the elution reaction plate on a magnet.

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

6.34. Using P200 multichannel set to 50µL volume, transfer all of the supernatant to a new 96 well plate. Make sure not to transfer any beads.

6.35. Cover the IP plate and keep at RT and proceed to the next step.

7. DNA Purification

7.1. Retrieve INPUT sample tube from 4°C storage.

7.2. Spin down at 200g for 1min at **ROOM TEMPERATURE** and place at RT.

7.3. Add 18µL of Qiagen EB buffer to each Input sample and mix slowly by pipetting up and down 10x. Final volume should be 30µL, to match the volume of the IPs.

7.4. Using a single channel P200 pipette, transfer the INPUT(s) to the IP plate as according to the plate layout worksheet. Cover the plate.

7.5. On ice, prepare the DNA Purification Master Mix. Scale up as calculated on worksheet.

Reagent	Volume (µL)/Rxn
Qiagen EB	31
Qiagen Buffer G2	8
Qiagen Protease	1.75
Total	40.75

Table 5: DNA Purification Master Mix

7.6. On ice, aliquot 40µL of DNA Purification Master Mix per row of samples plus 5µL dead volume into a new 96 well plate.

7.7. Using a P200 multichannel, add 40µL of DNA Purification Master Mix to IP plate (including INPUT wells). Mix slowing up and down 15x. Change tips after each row.

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 23 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

- 7.8. Seal the plate with Bio-Rad domed cap strips and incubate in a thermomixer at 50°C for 30 mins with a speed of 600 rpm (volume should be 70µL).
- 7.9. In preparation for Step 8, retrieve SeraMag beads and leave at RT for at least 30min in the dark. In addition, retrieve a new aliquot of 70% EtOH and keep at RT.
- 7.10. After incubation, spin plate at **RT** at 200g for 1 min.

8. Manual SeraMag Bead Clean Up of IPs

- 8.1. Retrieve SeraMag beads (30% PEG) from 4°C fridge and incubate at RT for at least 30 min before use. It is important that the beads are at room temperature as a cooler solution will yield lower product. Mix the beads well immediately before use. Due to the high percentage of PEG, the solution is quite viscous. Vortexing the tube is required to reach a homogenous suspension.
- 8.2. Retrieve 5mL of 70% EtOH and place at RT. Pour in a plate reservoir and cover for later use.
- 8.3. In a new plate, add 35µL of Qiagen EB buffer plus 5µL of dead volume per sample well. Cover and leave at RT for later use.
- 8.4. To the IP plate samples, add equal volume of beads (70µL). As long as the bead volume is equal or just a little more is sufficient.
- 8.5. Mix by pipetting up and down gently 10x.
- 8.6. Cover the plate and incubate at RT for 10-12min.
- 8.7. Put the plate on the magnet and let the solution become clear (approximately 10min)
- 8.8. Remove and discard the supernatant. Leave plate on magnet.
- 8.9. Add 150µL of 70% EtOH to each sample well.
- 8.10. Do not mix the beads and 70% EtOH as this is only a wash. Remove the supernatant.

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 24 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

- 8.11. Repeat 70% EtOH wash.
- 8.12. Take the plate off of the magnet.
- 8.13. Allow the beads to dry which takes a few minutes. (If the beads pellet appears cracked, then they are too dry).
- 8.14. Add 35µL of Qiagen EB buffer to the sample wells. Thoroughly mix by pipetting up and down at least 10x or until the mixture appears homogenous. Cover.
- 8.15. Incubate the plate at RT for 3min. Label a new 96 well plate with plate ID, date and initials and cover.
- 8.16. Place plate back on magnet and let sit for 2-3min.
- 8.17. Carefully transfer the supernatant to the new 96 well plate pre-labeled in Step 8.15.
- 8.18. Seal the plate with a foil cover and store at 4°C overnight or at -20°C for long term storage.

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

Appendix A - Buffer Preparation

NOTE: All buffers have a 2 month expiry unless otherwise noted

Reagent	Stock Concentration	Volume	Final Concentration
IP Buffer			
Tris-HCl (pH 7.5)	1M	200µL	20mM
EDTA	0.5M	40µL	2mM
NaCl	5M	300µL	150mM
Triton X-100**	1%	1mL	0.1%
Deoxycholate*	Solution		0.1%
Ultrapure Water	N/A	Up to 10mL	
Protease Inhibitor	Refer to protocol-add at the time of use		

Lysis Buffer			
Triton X-100**	1%	1mL	0.1%
Deoxycholate*	Solution		0.1%
Ultrapure Water	N/A	Up to 10mL	
Protease Inhibitor	Refer to protocol-add at the time of use		

Low Salt Wash Buffer			
Tris-HCl (pH 8.0)	1M	200µL	20mM
EDTA	0.5M	40µL	2mM
NaCl	5M	300µL	150mM
Triton-X 100**	10%	1mL	1%
SDS***	10%	100µL	0.1%
Ultrapure Water	N/A	Up to 10mL	

High Salt Wash Buffer			
Tris-HCl (pH 8.0)	1M	200µL	20mM
EDTA	0.5M	40µL	2mM
NaCl	5M	1mL	500mM
Triton-X 100**	10%	1mL	1%
SDS***	10%	100µL	0.1%
Ultrapure Water	N/A	Up to 10mL	
Elution Buffer			
NaHCO ₃ ****	1M****	100µL	100mM
SDS	10%	100µL	1%
Ultrapure Water	N/A	800µL	
Make a fresh solution for every use			

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 26 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

Stock Reagents:

***10% Triton X-100:**

1. Add 9mL of ultrapure water to a 15mL falcon tube.
2. Cut the end of a P1000 tip and add 1mL of Triton X-100.
3. Vortex/shake the tube until Triton X-100 is completely dissolved.
4. Label and date. Store at 4°C.
5. Expiry date is 2 months.

****1% Triton X-100+1% Deoxycholate:**

1. Using analytical scale, weigh out 100mg of sodium deoxycholate into a 15mL falcon tube.
2. Add 1mL of 10% Triton X-100*.
3. Top up to 10mL with ultrapure water.
4. Vortex/shake vigorously until all components are dissolved.
5. Label and date. Store at 4°C.
6. Expiry date is 2 months.

*****10% SDS:**

1. Add 5mL of 20% SDS solution to a 15mL falcon tube.
2. Top up to 10mL with ultrapure water.
3. Mix gently to resuspend.
4. Label and date. Store at room temperature.
5. Expiry is 6 months.

******1M NaHCO₃:**

1. Using an analytical scale, weigh out 840mg of NaHCO₃ into a 15mL falcon tube.
2. Top up to a final volume of 10mL with ultrapure water.
3. Vortex the solution vigorously until NaHCO₃ is completely dissolved.
4. Label and date. Store at room temperature.
5. Expiry date is 1 month.

Note: All buffers, enzymes will need to be validated using the Native ChIP production pipeline up to library qPCR QC.

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 27 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

Appendix B-Dilution of Stock Antibodies

-if the concentration of the stock is too high to pipette a minimum of 2µL/IP then it requires a dilution.

-the buffer to use to dilute the stock is the same buffer used to make up the stock, which is provided on the information sheet provided by the supplier with the shipment of the antibody.

-for example, Figure 4 details the components of the storage buffer of the H3K4me1 antibody

H3K4me1 polyclonal antibody

<p>Cat. No. C15410194 (pAb-194-050)</p> <p>Type: Polyclonal ChIP-grade, ChIP-seq grade</p> <p>Source: Rabbit</p> <p>Lot #: A1862D</p> <p>Size: 50 µg/ 34 µl</p> <p style="border: 1px solid red; padding: 2px;">Concentration: 1.5 µg/µl</p> <p style="text-align: center;">↑</p> <p>Concentration of Stock Antibody</p>	<p>Specificity: Human, mouse, wide range: expected</p> <p style="border: 1px solid red; padding: 2px;">Purity: Affinity purified polyclonal antibody in PBS containing 0.05% azide and 0.05% ProClin 300.</p> <p>Storage: Store at -20°C; for long storage, store at -80°C. Avoid multiple freeze-thaw cycles.</p> <p>Precautions: This product is for research use only. Not for use in diagnostic or therapeutic procedures.</p>
--	--

← Storage buffer Details

Figure 4

- For this particular antibody the concentration is 1.5µg/uL.
- The amount of antibody needed per IP is 0.5µg, which is 0.33µL. This is too low to pipette.
- It is recommended to take an aliquot of the antibody and dilute to a concentration that is easily pipettable.
- Consult with supervisor if unsure.

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 28 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

Appendix C-Preparation of 30% PEG SeraMAG Beads

1. 30% PEG in 1M NaCL buffer preparation

- 1.1. Place a new, sterile 50mL falcon tube on the analytical scale plate (use a styrofoam rack) and zero it. Accurately weigh 15.0g of PEG-8000 directly into the 50 mL tube. Repeat when preparing multiple tubes.
- 1.2. Add 20mL of ultrapure water (DNase/RNase free distilled water). Close the tube and mix by inverting.
- 1.3. Add 10mL of 5M NaCL.
- 1.4. Add 500µL of 1M Tris-HCL (pH8.0).
- 1.5. Add 100µL of 0.5M EDTA.
- 1.6. Mix by inverting until PEG goes into the solution and all of the components are uniformly dispersed. The solution should become completely clear. Use a nutating mixer if preparing more than one tube at a time (recommended).
- 1.7. Label the tube with a permanent marker and keep the mixture at room temperature while washing the beads.

2. Bead Washes

- 2.1. Remove SeraMag beads from its 4°C storage, mix very well to make sure that the beads are evenly dispersed, aliquot 1mL into a new 1.5mL non-stick tube.
- 2.2. Place the tube in the DynaMag magnet stand and incubate for 2min. Make sure all the beads are drawn to the sides and the solution is completely clear.
- 2.3. Remove the supernatant and discard. Be very careful not to remove any beads. Take the tube off the magnet stand.
- 2.4. Add 1mL of TE buffer to the beads, close the tube and mix by gentle, repeated pulse-vortexing.
- 2.5. Place the tube back on the magnet and incubate for 2min. Make sure that all of the beads separate to the side and the solution becomes completely clear.

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 29 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

- 2.6. Remove and discard the supernatant. Be careful not to remove any beads.
- 2.7. Take the tube off the magnet and repeat TE buffer wash two more times (for a total of 3 washes).
- 2.8. After the last wash resuspend the beads in 1mL of TE buffer and keep it at room temperature.

3. SeraMag Bead Stock Preparation

- 3.1. Aliquot 1mL of the washed SeraMag beads to the 30% PEG, 1M NaCl solution. Resuspend the beads by inversion.
- 3.2. Top up the solution to 50 mL with ultrapure water, close the lid and mix by inverting and *gentle*, repeated, pulse vortexing until uniformly brown.
- 3.3. When preparing multiple 50mL conical tubes of beads at a time, mix all vials together in a pre-sterilized glass bottle to create one batch solution. This reduces the number of QC tests.
- 3.4. Label the tube (or the glass bottle) as “Do not use. QC in progress.”
- 3.5. Remove an aliquot of the freshly prepared SeraMag bead solution for QC testing. Store the rest at 4°C, **protected from light**.

4. Quality Control Testing

- 4.1. Test every new batch of beads prepared by performing Native ChIP using 100,000 HL60 cells. IP whole histone mod panel.
- 4.2. Perform library construction including Quant-iT QC for product amount and Caliper for size distribution.
- 4.3. Perform qPCR using histone modification target specific primers and calculate fold enrichment.
- 4.4. Compile all results and confer with supervisor.
- 4.5. If the bead batch passes all QC gates proceed to Step.5. If batch production was unsuccessful prepare a fresh batch.

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 30 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

5. Aliquot Beads

- 5.1. After the new batch passes QC, retrieve the newly prepared SeraMag bead working stock from its designated 4°C storage. Make sure the stock fully reaches room temperature (for larger volumes this can take 1h or more).
- 5.2. Mix the SeraMag bead stock very well by inverting and gentle, repeated, pulse-vortexing. Make sure that the beads are uniformly dispersed in the supernatant and that the solution is uniform in color.
- 5.3. Aliquot 1mL of beads to 1.5mL non-stick prelabelled tubes.
- 5.4. Store aliquots at 4°C, **protected from light**. Make sure to label and date each new batch on the tubes and on the storage box.
- 5.5. Bead solution stored at 4°C and protected from light should be stable for at least 6 months. After 6 months re-test monthly for performance.

Native ChIP Using 100,000 Cells	
Document #: LIBPR.0138	Supersedes: Version 1
Version: 2	Page 31 of 31

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

Appendix D:

LIMS:

1. Native ChIP
2. Native ChIP qPCR (to be done after library construction and qPCR QC)

OFFICIAL