

Operation and Maintenance of the Agilent 2100 Bioanalyzer for RNA Samples	
Saved as: LIBPR.0018	Supersedes: Version 8
Version: 9	Page 1 of 19

## Non Controlled Version

*\*Note: Controlled Versions of this document are subjected to change without notice*

# Operation and Maintenance of the Agilent 2100 Bioanalyzer for RNA Samples

## I. Purpose

To provide specific guidelines for operating and maintaining the Agilent 2100 Bioanalyzer. This is achieved using reagents and instructions provided by Agilent to determine the quality and quantity of RNA samples.

## II. Scope

All procedures are applicable to the BCGSC FG-Library Construction and TechD groups.

## III. Policy

All production procedures shall be documented and controlled by approved systems.

## IV. Responsibility

It is the responsibility of all personnel performing this procedure to follow the current protocol. It is the responsibility of the Group Leader to ensure personnel are trained in all aspects of this protocol. It is the responsibility of Quality Assurance Management to audit this procedure for compliance and maintain control of this procedure.

## V. References

Reference Title	Reference Number
Agilent RNA 6000 Nano Assay Quick Start Guide	G2938-90030
Agilent RNA 6000 Pico Assay Quick Start Guide	G2941-90161

## VI. Related Documents

Document Title	Document Number
N/A	N/A

## VII. Safety

All Laboratory Safety procedures will be complied with during this procedure. The required personal protective equipment includes a laboratory coat and gloves. See the material safety data sheet (MSDS) for additional information.

Operation and Maintenance of the Agilent 2100 Bioanalyzer for RNA Samples	
Saved as: LIBPR.0018	Supersedes: Version 8
Version: 9	Page 2 of 19

## Non Controlled Version

\*Note: Controlled Versions of this document are subjected to change without notice

### VIII. Materials and Equipment

Name	Supplier	Number: #	Model or Catalogue #
Small Autoclave waste bags 10”X15”	Fisher Scientific	01-826-4	✓
RNase Zap	Ambion	9780	✓
Ice bucket – Blue	Fisher	11-676-36	✓
wet ice	In house	N/A	N/A
Gilson P2 pipetman	Mandel	GF-44801	✓
Gilson P10 pipetman	Mandel	GF-44802	✓
Gilson P1000 pipetman	Mandel	GF-23602	✓
Diamond Filter Tips 10ul	Mandel	GF-F171203	✓
Diamond Filter Tips 1000 µl	Mandel	GF-F171703	✓
Galaxy mini-centrifuge	VWR	37000-700	✓
VX-100 Vortex Mixer	Rose Scientific	S-0100	✓
Large Kimwipes	Fisher	06-666-117	✓
Black ink permanent marker pen	VWR	52877-310	✓
Eppendorf BenchTop Refrigerated Centrifuge 5810R	Eppendorf	5810 R	✓
Agilent RNA 6000 Pico Kit Series II	Agilent	5067-1513	✓
Agilent RNA 6000 Pico Kit Reagents Series II	Agilent	5067-1514	✓
Agilent RNA 6000 Nano Kit Series II	Agilent	5067-1511	✓
Agilent RNA 6000 Nano Reagents Series II	Agilent	5067-1512	✓
Chip Priming Station	Agilent	5065-4401	✓
IKA Vortex Mixer	Agilent	MS2-S8	✓
Agilent 2100 Bioanalyzer	Agilent	2100	✓
Computer			✓
Heat Block ISOTEMP 125D	Fisher	125D	✓
RNase Zap	Ambion	AM9780	✓
DEPC water	Ambion	AM9920	✓
1.5ml Eppendorf Tubes	Ambion	AM12400	

### IX. Procedure

**Note: All Chip priming, running of the Agilent 2100 Bioanalyzer is to be completed on the 6<sup>th</sup> Floor. (Figure 1)**

Operation and Maintenance of the Agilent 2100 Bioanalyzer for RNA Samples	
Saved as: LIBPR.0018	Supersedes: Version 8
Version: 9	Page 3 of 19

## Non Controlled Version

\*Note: Controlled Versions of this document are subjected to change without notice



**Figure 1. Agilent 2100 Bioanalyzer System**

### 1. Determine Which Assay Is Needed Based On Sample Type And Expected Amount Of Material

- 1.1. RNA samples can either be analyzed on RNA 6000 Nano or RNA 6000 Pico chip.
- 1.2. With the expected amount of RNA, select the appropriate assay based on the following table:

Specification	Total RNA Nano Assay	Total RNA Pico Assay
Quantitative Range	25-500ng/μl	n/a
Qualitative Range	5-500ng/μl	50-5000pg/μl

**\*Note:** mRNA assay is only for qualitative QC, not quantitative. Use 1/5<sup>th</sup> dilution of the RNA Nano marker for the mRNA Agilent assay in step 7.2 if less than 10μg of total RNA was used for mRNA purification. If using the diluted RNA Nano marker, a threshold of 6% of rRNA contamination is acceptable, discuss with supervisor if rRNA contamination % is above 6%. mRNA samples may also be run using the “Total Eukaryotic RNA assay” setting and regular RNA Nano Marker without dilution; however, the % of rRNA contamination calculation is not a feature of this assay.

### 2. Replacing the Chip Priming Station Syringe

Note: This step is done as a weekly maintenance procedure by Library Core. It does not need to be performed each time a chip is being primed.

- 2.1. Put on a clean pair of gloves

Operation and Maintenance of the Agilent 2100 Bioanalyzer for RNA Samples	
Saved as: LIBPR.0018	Supersedes: Version 8
Version: 9	Page 4 of 19

## *Non Controlled Version*

*\*Note: Controlled Versions of this document are subjected to change without notice*

- 2.2. Unscrew the old syringes from the lid of the Chip Priming Station (see Figure 2)
- 2.3. Remove the plastic cap of the new syringe and insert it into the clip
- 2.4. Slide it into the hole of the luer lock adapter and screw it tightly to the Chip Priming Station




Figure 2 Chip Priming Station

### **3. Filtered Gel Matrix Preparation**

- 3.1. Remove the reagents for the desired kit from the 4°C fridge (EQU2014) and allow them to equilibrate to room temperature for 30 minutes before use. Put reagents in pre-made box that shelters dye and gel-dye mix from light. RNA ladders for both Pico and Nano kits are kept in a box in the -80°C freezer (EQU275). If unsure of where reagents are kept, please check with supervisor.
- 3.2. Pipette 550ul of RNA Nano or Pico gel matrix into a spin filter.
- 3.3. Centrifuge at 4000rpm for 10 minutes at room temperature.
- 3.4. Remove the filter and label the tube as “Filtered Gel Matrix” for either Nano or Pico Assay, with the date and your initials
- 3.5. Use the filtered gel within 4 weeks. Discard afterwards



## Non Controlled Version

*\*Note: Controlled Versions of this document are subjected to change without notice*

### 4. Gel Dye Mix Preparation for Nano and Pico chips

- 4.1. Spin the Gel Matrix briefly to collect at the bottom of the tube.
- 4.2. Transfer 65  $\mu\text{L}$  of the filtered gel matrix into a fresh RNase-free 1.5 mL tube. Label the tube as either "Nano Gel Dye" or "Pico Gel Dye" with the date and initials.
- 4.3. Vortex the dye concentrate for 10 seconds at high speed. Briefly spin down to collect solution at the bottom of the tube.
- 4.4. Add 1  $\mu\text{L}$  of the dye concentrate to the 65  $\mu\text{L}$  of filtered gel matrix. Return to dye concentrate to the box, shielded from light.
- 4.5. Vortex the solution at high speed for 10 seconds. Spin at 13,000 rpm for 10 minutes at room temperature. Keep tube upright and in the dark.

### 5. Setting Up and Running the Chip

- 5.1. Put on a clean pair of gloves.
- 5.2. Aliquot and dilute samples to the appropriate range while reagents are equilibrating. If an aliquot is to be made, 1.5  $\mu\text{L}$  of sample is required. Keep all samples on ice. Arrange samples by approximate lowest to highest concentration. Arrange in the order as they will be loaded onto the chip to prevent mix-up. **To remove any secondary structure place samples and the RNA Nano ladder onto 70°C heat block for 2 minutes and then snap cool samples by returning tubes to ice for another 2 minutes followed by a pulse spin.**

**\*Note:** It is NOT necessary to heat the Pico Ladder as it has already been denatured.

- 5.3. Ensure that the Chip Priming Station is set up correctly for type of chip being run. The base plate should be in position "C".
- 5.4. If running a **Nano** chip assay, clean the electrode before and after as follows:
  - 5.4.1. Slowly fill one of the wells of the electrode cleaner chip with 350  $\mu\text{L}$  of fresh RNase Zap

Operation and Maintenance of the Agilent 2100 Bioanalyzer for RNA Samples	
Saved as: LIBPR.0018	Supersedes: Version 8
Version: 9	Page 6 of 19

## *Non Controlled Version*

*\*Note: Controlled Versions of this document are subjected to change without notice*

- 5.4.2. Never fill too much fluid in the electrode cleaner. This could cause liquid spill which may cause current leaks between electrodes.
  - 5.4.3. Open the lid and place electrode cleaner chip in the Agilent 2100 Bioanalyzer.
  - 5.4.4. Close the lid and leave it closed for 1 minute.
  - 5.4.5. Open the lid and remove the electrode cleaner chip.
  - 5.4.6. Slowly fill one of the wells of another electrode cleaner chip with 350 µL RNase Free water.
  - 5.4.7. Insert the RNase free water cleaner chip into the Agilent 2100 Bioanalyzer and close the lid. Leave it closed for 30 seconds.
  - 5.4.8. Open the lid and remove the electrode cleaner chip. Leave the lid open for 30 seconds to allow the water on the electrodes to evaporate.
  - 5.4.9. Make sure there is a date on the cleaner chips and keep for the end of protocol cleaning. Following the end of protocol cleaning, tap out the liquid of each cleaner chip and leave upside down on a paper towel adjacent to the Bioanalyzer
- 5.5. If running a **Pico** chip assay, clean the electrode before and after as follows.
- 5.5.1. Slowly fill one of the wells of the electrode cleaner chip with 350µL of fresh RNase free water.
  - 5.5.2. Never fill too much fluid in the electrode cleaner. This could cause liquid spill which might cause current leaks between the electrodes.
  - 5.5.3. Open the lid and place the electrode cleaner chip in the Agilent 2100 Bioanalyzer.
  - 5.5.4. Close the lid and leave it closed for 5 minutes.
  - 5.5.5. Open the lid and remove the electrode cleaner chip. Leave the lid open and wait for 30 seconds for the water on the electrodes to evaporate.
  - 5.5.6. Make sure there is a date on the cleaner chip and keep for future use.
- 5.6. Ensure that the Chip Priming Station is set up correctly for type of chip being run. The base plate should be in position “C” and the clip should be in the topmost position.
- 5.7. Turn on the Agilent 2100 Bioanalyzer by flipping the switch at the rear of the instrument. A green light will come on.
- 5.7.1. Login to the computer supporting the Agilent 2100 Bioanalyzer and double click on the “2100 Expert” icon.
  - 5.7.2. The screen of the software appears in the *Instrument* context.
  - 5.7.3. Click on the *Assay Selection* icon and choose the correct assay to be run. Enter the number of samples to be run in the appropriate box.
  - 5.7.4. Under “Destination” click on the circle adjacent to “Custom” and save the file
    - 5.7.4.1. For Nano under *R:\ Library Core\QC\Agilent\RNA Nano*
    - 5.7.4.2. For Pico under *R:\ Library Core\QC\Agilent\RNA*

Operation and Maintenance of the Agilent 2100 Bioanalyzer for RNA Samples	
Saved as: LIBPR.0018	Supersedes: Version 8
Version: 9	Page 7 of 19

## *Non Controlled Version*

*\*Note: Controlled Versions of this document are subjected to change without notice*

- 5.7.4.3. For mRNA under *R:\ Library Core\QC\Agilent\mRNA*
- 5.7.5. Identify the chip with the next assay# by entering this new *assay#* in the file prefix section (refer to whiteboard to naming scheme: **RNA#** (Assay) for Nano; **Pico#** (Assay) for Pico; **mRNA#** (Assay) for mRNA.
- 5.7.6. Change the Assay number on the whiteboard to the Assay number currently being run.

## 6. Chip Priming

- 6.1. Set a timer for 30 seconds.
- 6.2. Remove a new RNA Nano or Pico Chip from its sealed bag.
- 6.3. When pipetting the gel dye mix, be careful not to draw up any particles that may sit at the bottom of the tube. Pipette 9.0  $\mu$ L of the gel dye mix into the bottom of the well marked G (with black circle). Take care not to introduce bubbles when dispensing the gel dye mix. If bubbles are present, very gently tap the chip on the bench top or remove the gel dye mix from the well and repipette 9.0  $\mu$ L. The chip will not run properly if air gaps are pushed into the capillaries.




Figure 3. Proper placement of pipette tip in chip well.

- 6.4. Place the chip onto the Chip Priming station.
- 6.5. Raise the syringe plunger to the 1 mL mark. Ensure the lever clip is in the top position.
- 6.6. Press the plunger until it is held by the clip.
- 6.7. Start the timer, wait exactly 30 seconds then release the clip. Allow the plunger to rise by itself.
- 6.8. Wait for 5 seconds before slowly pulling the plunger back to the 1 mL mark.
- 6.9. Remove the chip from the Chip Priming station.

Operation and Maintenance of the Agilent 2100 Bioanalyzer for RNA Samples	
Saved as: LIBPR.0018	Supersedes: Version 8
Version: 9	Page 8 of 19

## *Non Controlled Version*

*\*Note: Controlled Versions of this document are subjected to change without notice*

### **7. Loading the Remaining Gel Dye Mix and Nano/Pico Marker**

- 7.1. Pipette 9.0 µL of gel dye mix into the 2 wells marked G.
- 7.2. Vortex the RNA Nano or Pico Marker for 3 seconds at high speed then quick spin. If running mRNA assay to check mRNA purified from less than 10µg of total RNA, dilute the RNA Nano marker by 1/5<sup>th</sup>. Pipette 5.0µL of the RNA Nano or diluted RNA Nano or Pico Marker, depending on which assay it is, into the ladder well and then into all sample wells. Do not leave any empty wells or the chip will not run.
- 7.3. For Pico chip assay, pipette 9ul of the RNA 6000 Pico conditioning solution into the well marked CS.

### **8. Loading the Ladder and samples**

***Note: Ensure the samples and RNA Nano ladder have been denatured for 2 minutes at 70 degrees and then snap chilled for 2 minutes before addition onto the chip (See 3.3)***

- 8.1. Pipette 1.0 µL of denatured RNA Nano or Pico Ladder into the ladder well.
- 8.2. Pipette 1.0 µL of denatured sample into each of the sample wells.
- 8.3. Pipette 1.0 µL of DEPC-water or RNA Marker into each unused sample well.
- 8.4. Put the chip in the Agilent Vortex Mixer adapter and secure with labeling tape if necessary. Set speed at **2200 rpm** and vortex chip for one minute. ***A prepared chip must be loaded within 5 minutes.***
- 8.5. Open the lid of the Bioanalyzer and carefully place the chip into the receptacle. The chip fits only one way.
- 8.6. Carefully close the lid. The electrodes fit into the well of the chip. The 2100 Expert software should display the appropriate assay chip on the screen.
- 8.7. Ensure that the appropriate assay has been selected.
- 8.8. Enter the sample information, including the Library ID or Sample ID, dilution information, volume of RNA sample and Initials in the chip summary table.



Operation and Maintenance of the Agilent 2100 Bioanalyzer for RNA Samples	
Saved as: LIBPR.0018	Supersedes: Version 8
Version: 9	Page 9 of 19

## Non Controlled Version

\*Note: Controlled Versions of this document are subjected to change without notice

8.9. Click the *Start* button in the upper right window. The incoming raw signals are displayed in the *Instrument* context.




Figure 4. Sample of a Nano Ladder.




Figure 5. Sample of a Pico Ladder

## 9. Shutdown and cleaning



Operation and Maintenance of the Agilent 2100 Bioanalyzer for RNA Samples	
Saved as: LIBPR.0018	Supersedes: Version 8
Version: 9	Page 10 of 19

## *Non Controlled Version*

*\*Note: Controlled Versions of this document are subjected to change without notice*

- 9.1. Immediately after the chip is finished, remove the chip from the receptacle of the Bioanalyzer and dispose of it in the waste container adjacent the instrument.
- 9.2. Refer to Steps 5.4 and 5.5 for the appropriate cleaning procedure.
- 9.3. Close the software and Log off of the computer. Turn off the monitor as well. Switch off the Agilent 2100 BioAnalyzer.

OFFICIAL

Operation and Maintenance of the Agilent 2100 Bioanalyzer for RNA Samples	
Saved as: LIBPR.0018	Supersedes: Version 8
Version: 9	Page 11 of 19


## Non Controlled Version

\*Note: Controlled Versions of this document are subjected to change without notice

### Appendix A

#### 1. Creating Bioanalyzer Runs for you Samples

1.1. Under the Lib\_Construction tab, scan in your samples and the Agilent equ# that you used (equ283 or equ2216).



1.2. A list will be generated of all the samples that you have scanned in.

Create Bioanalyzer Run for the following plates

Plate ID	Library	Equipment Name
<a href="#">678752</a>	<a href="#">A29912:CD_1428</a>	BA2100-7
<a href="#">678753</a>	<a href="#">A29913:CD_1429</a>	BA2100-7
<a href="#">678754</a>	<a href="#">A29914:CD_1430</a>	BA2100-7
<a href="#">678755</a>	<a href="#">A29915:CD_1431</a>	BA2100-7
<a href="#">678756</a>	<a href="#">A29916:CD_1432</a>	BA2100-7
<a href="#">678757</a>	<a href="#">A29917:CD_1433</a>	BA2100-7
<a href="#">678758</a>	<a href="#">A29918:CD_1434</a>	BA2100-7
<a href="#">678759</a>	<a href="#">A29919:CD_1435</a>	BA2100-7
<a href="#">678760</a>	<a href="#">A29920:CD_1436</a>	BA2100-7
<a href="#">678761</a>	<a href="#">A29921:CD_1437</a>	BA2100-7
<a href="#">678762</a>	<a href="#">A29922:CD_1438</a>	BA2100-7

11 Records

Sample positions	
Plate	Position(s)
678752	<input type="text"/>
678753	<input type="text"/>

Operation and Maintenance of the Agilent 2100 Bioanalyzer for RNA Samples	
Saved as: LIBPR.0018	Supersedes: Version 8
Version: 9	Page 12 of 19

## Non Controlled Version

\*Note: Controlled Versions of this document are subjected to change without notice

- 1.3. In the “Sample Positions” table, fill in the position each sample loaded on the Agilent chip. Even though it is a 4x4 well chip, only the sample wells will be ‘numbered’. So in row 1 of the chip, you have positions 1-3. Row 2 has positons 4-6 and so on. For the mRNA and RNAnano assays, there is a maximum of 12 postions. For the RNApico assay, there is a maximum of 11 postions.

Sample positions	
Plate	Position(s)
678752	<input type="text"/>
678753	2 <input type="text"/>
678754	3 <input type="text"/>
678755	4 <input type="text"/>
678756	5 <input type="text"/>
678757	6 <input type="text"/>
678758	7 <input type="text"/>
678759	8 <input type="text"/>
678760	9 <input type="text"/>
678761	10 <input type="text"/>
678762	11 <input type="text"/>

Enter position(s) of this sample on the chip. e.g. 1,2,3 or 1-3

- 1.4. Enter in the Chip ID. If you accidentally enter in a Chip ID that was previously used for this assay, the next available Chip ID number will be assigned to your samples. Select, “Create Bioanalyzer Run”.

Operation and Maintenance of the Agilent 2100 Bioanalyzer for RNA Samples	
Saved as: LIBPR.0018	Supersedes: Version 8
Version: 9	Page 13 of 19

## Non Controlled Version

\*Note: Controlled Versions of this document are subjected to change without notice

678760	<input type="text" value="9"/>
678761	<input type="text" value="10"/>
678762	<input type="text" value="11"/>

Input Parameters

Chip ID:

Comments:

Create Bioanalyzer Run

1.5. One run is created for each sample you scanned in.

Creating Bioanalyzer Run(s)

Recorded Storage Location for Plates; reset to 'In Use'

Recorded relocation

Moved 1 Plate(s) -> [Rac10\\_In Use](#)

Created Run(s): [133314](#), [133315](#), [133316](#), [133317](#), [133318](#), [133319](#), [133320](#), [133321](#), [133322](#), [133323](#), [133324](#) Please go to 'Bioanalyzer Runs' tab to monitor the progress and to view the run results.

1.6. To monitor the progress and to view the run results, select the “Bioanalyzer Runs” tab.

Lib\_Construction Receiving

Scan Grab Plate Set Error Search DB

[Site Map](#) [Lab](#) [Admin](#) [Database](#) [Shipments](#) [Views](#) [Summaries](#) [Refresh Settings](#)

---

Creating Bioanalyzer Run(s)

Recorded Storage Location for Plates; reset to 'In Use'

Recorded relocation

Moved 1 Plate(s) -> [Rac10\\_In Use](#)

Created Run(s): [133295](#), [133296](#), [133297](#), [133298](#), [133299](#), [133300](#), [133301](#), [133302](#), [133303](#), [133304](#), [133305](#) Please go to 'Bioanalyzer Runs' tab to monitor the progress and to view the run results.

Current Samples: [1 Tube](#)

Lab
Bioanalyzer Runs
Admin

1.7. Your samples will show up under “Pending Runs”.

Operation and Maintenance of the Agilent 2100 Bioanalyzer for RNA Samples	
Saved as: LIBPR.0018	Supersedes: Version 8
Version: 9	Page 14 of 19

## Non Controlled Version

\*Note: Controlled Versions of this document are subjected to change without notice

Lab
Bioanalyzer Runs
Admin

---

**Bioanalyzer Runs**

No Waiting Analysis runs found;

---

Pending Runs														
Position	Plate	Library	Run_ID	Run_Status	Run_DateTime	Run_Finished	Result_File	Chip_ID	Assay_Type	Equipment	Run_Test_Status	Run_Validation	Run_QC_Status	Run_Comments
1	678752	A29912:CD_1428	133314	In Process	Apr-19-2013 15:28:48	-	-	RNA#1550	-	BA2100-7 - EQU2216	Production	Pending	N/A	-
2	678753	A29913:CD_1429	133315	In Process	Apr-19-2013 15:28:48	-	-	RNA#1550	-	BA2100-7 - EQU2216	Production	Pending	N/A	-
3	678754	A29914:CD_1430	133316	In Process	Apr-19-2013 15:28:48	-	-	RNA#1550	-	BA2100-7 - EQU2216	Production	Pending	N/A	-
4	678755	A29915:CD_1431	133317	In Process	Apr-19-2013 15:28:48	-	-	RNA#1550	-	BA2100-7 - EQU2216	Production	Pending	N/A	-
5	678756	A29916:CD_1432	133318	In Process	Apr-19-2013 15:28:48	-	-	RNA#1550	-	BA2100-7 - EQU2216	Production	Pending	N/A	-
6	678757	A29917:CD_1433	133319	In Process	Apr-19-2013 15:28:48	-	-	RNA#1550	-	BA2100-7 - EQU2216	Production	Pending	N/A	-
7	678758	A29918:CD_1434	133320	In Process	Apr-19-2013 15:28:48	-	-	RNA#1550	-	BA2100-7 - EQU2216	Production	Pending	N/A	-
8	678759	A29919:CD_1435	133321	In Process	Apr-19-2013 15:28:48	-	-	RNA#1550	-	BA2100-7 - EQU2216	Production	Pending	N/A	-
9	678760	A29920:CD_1436	133322	In Process	Apr-19-2013 15:28:48	-	-	RNA#1550	-	BA2100-7 - EQU2216	Production	Pending	N/A	-
10	678761	A29921:CD_1437	133323	In Process	Apr-19-2013 15:28:48	-	-	RNA#1550	-	BA2100-7 - EQU2216	Production	Pending	N/A	-
11	678762	A29922:CD_1438	133324	In Process	Apr-19-2013 15:28:48	-	-	RNA#1550	-	BA2100-7 - EQU2216	Production	Pending	N/A	-
11 Records														

Analyzed Runs

1.8. Once the Runs have been validated by LIMS, your samples will show up in the “Waiting Analysis Runs” section. LIMS will edit the Chip ID number if incorrect. It will also fill in when the run was completed, where the results are stored, and the assay type. This takes about 10 minutes.

## Non Controlled Version

\*Note: Controlled Versions of this document are subjected to change without notice

Lab
Bioanalyzer Runs
Admin

---

**Bioanalyzer Runs**

Waiting Analysis Runs

Toggle	Position	Plate	Library	Run_ID	Run_Status	Run_DateTime	Run_Finished	Result_File	Chip_ID	Assay_Type	Equipment	Run_Test_Status	Run_Validation	Run_QC_Status	Run_Comments
<input type="checkbox"/>	1	678752	<a href="#">A29912:CD_1428</a>	<a href="#">133314</a>	Data Acquired	Apr-18-2013 00:00:00	2013-04-19 15:41:22	\\isaac\binstruments\Bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	<a href="#">BA2100-4-EQU283</a>	Production	Pending	N/A	-
<input type="checkbox"/>	2	678753	<a href="#">A29913:CD_1429</a>	<a href="#">133315</a>	Data Acquired	Apr-18-2013 00:00:00	2013-04-19 15:41:22	\\isaac\binstruments\Bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	<a href="#">BA2100-4-EQU283</a>	Production	Pending	N/A	-
<input type="checkbox"/>	3	678754	<a href="#">A29914:CD_1430</a>	<a href="#">133316</a>	Data Acquired	Apr-18-2013 00:00:00	2013-04-19 15:41:22	\\isaac\binstruments\Bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	<a href="#">BA2100-4-EQU283</a>	Production	Pending	N/A	-
<input type="checkbox"/>	4	678755	<a href="#">A29915:CD_1431</a>	<a href="#">133317</a>	Data Acquired	Apr-18-2013 00:00:00	2013-04-19 15:41:22	\\isaac\binstruments\Bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	<a href="#">BA2100-4-EQU283</a>	Production	Pending	N/A	-
<input type="checkbox"/>	5	678756	<a href="#">A29916:CD_1432</a>	<a href="#">133318</a>	Data Acquired	Apr-18-2013 00:00:00	2013-04-19 15:41:22	\\isaac\binstruments\Bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	<a href="#">BA2100-4-EQU283</a>	Production	Pending	N/A	-
<input type="checkbox"/>	6	678757	<a href="#">A29917:CD_1433</a>	<a href="#">133319</a>	Data Acquired	Apr-18-2013 00:00:00	2013-04-19 15:41:22	\\isaac\binstruments\Bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	<a href="#">BA2100-4-EQU283</a>	Production	Pending	N/A	-
<input type="checkbox"/>	7	678758	<a href="#">A29918:CD_1434</a>	<a href="#">133320</a>	Data Acquired	Apr-18-2013 00:00:00	2013-04-19 15:41:22	\\isaac\binstruments\Bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	<a href="#">BA2100-4-EQU283</a>	Production	Pending	N/A	-
<input type="checkbox"/>	8	678759	<a href="#">A29919:CD_1435</a>	<a href="#">133321</a>	Data Acquired	Apr-18-2013 00:00:00	2013-04-19 15:41:22	\\isaac\binstruments\Bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	<a href="#">BA2100-4-EQU283</a>	Production	Pending	N/A	-
<input type="checkbox"/>	9	678760	<a href="#">A29920:CD_1436</a>	<a href="#">133322</a>	Data Acquired	Apr-18-2013 00:00:00	2013-04-19 15:41:22	\\isaac\binstruments\Bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	<a href="#">BA2100-4-EQU283</a>	Production	Pending	N/A	-
<input type="checkbox"/>	10	678761	<a href="#">A29921:CD_1437</a>	<a href="#">133323</a>	Data Acquired	Apr-18-2013 00:00:00	2013-04-19 15:41:22	\\isaac\binstruments\Bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	<a href="#">BA2100-4-EQU283</a>	Production	Pending	N/A	-
<input type="checkbox"/>	11	678762	<a href="#">A29922:CD_1438</a>	<a href="#">133324</a>	Data Acquired	Apr-18-2013 00:00:00	2013-04-19 15:41:22	\\isaac\binstruments\Bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	<a href="#">BA2100-4-EQU283</a>	Production	Pending	N/A	-

11 Records

Add Analysis Data



1.9. Select the samples you want to analyze and select “Add Analysis Data”.

Operation and Maintenance of the Agilent 2100 Bioanalyzer for RNA Samples	
Saved as: LIBPR.0018	Supersedes: Version 8
Version: 9	Page 16 of 19

## Non Controlled Version

\*Note: Controlled Versions of this document are subjected to change without notice

Lab
Bioanalyzer Runs
Admin

Bioanalyzer Runs															
Waiting Analysis Runs															
Toggle	Position	Plate	Library	Run_ID	Run_Status	Run_DateTime	Run_Finished	Result_File	Chip_ID	Assay_Type	Equipment	Run_Test_Status	Run_Validation	Run_QC_Status	Run_Comments
<input checked="" type="checkbox"/>	1	678752	<a href="#">A29912:CD_1428</a>	<a href="#">133314</a>	Data Acquired	Apr-18-2013 00:00:00	2013-04-19 15:41:22	\\isaac\labinstruments\Bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	BA2100-4-EQU283	Production	Pending	N/A	-
<input checked="" type="checkbox"/>	2	678753	<a href="#">A29913:CD_1429</a>	<a href="#">133315</a>	Data Acquired	Apr-18-2013 00:00:00	2013-04-19 15:41:22	\\isaac\labinstruments\Bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	BA2100-4-EQU283	Production	Pending	N/A	-
<input checked="" type="checkbox"/>	3	678754	<a href="#">A29914:CD_1430</a>	<a href="#">133316</a>	Data Acquired	Apr-18-2013 00:00:00	2013-04-19 15:41:22	\\isaac\labinstruments\Bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	BA2100-4-EQU283	Production	Pending	N/A	-
<input checked="" type="checkbox"/>	4	678755	<a href="#">A29915:CD_1431</a>	<a href="#">133317</a>	Data Acquired	Apr-18-2013 00:00:00	2013-04-19 15:41:22	\\isaac\labinstruments\Bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	BA2100-4-EQU283	Production	Pending	N/A	-
<input checked="" type="checkbox"/>	5	678756	<a href="#">A29916:CD_1432</a>	<a href="#">133318</a>	Data Acquired	Apr-18-2013 00:00:00	2013-04-19 15:41:22	\\isaac\labinstruments\Bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	BA2100-4-EQU283	Production	Pending	N/A	-
<input checked="" type="checkbox"/>	6	678757	<a href="#">A29917:CD_1433</a>	<a href="#">133319</a>	Data Acquired	Apr-18-2013 00:00:00	2013-04-19 15:41:22	\\isaac\labinstruments\Bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	BA2100-4-EQU283	Production	Pending	N/A	-
<input checked="" type="checkbox"/>	7	678758	<a href="#">A29918:CD_1434</a>	<a href="#">133320</a>	Data Acquired	Apr-18-2013 00:00:00	2013-04-19 15:41:22	\\isaac\labinstruments\Bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	BA2100-4-EQU283	Production	Pending	N/A	-
<input checked="" type="checkbox"/>	8	678759	<a href="#">A29919:CD_1435</a>	<a href="#">133321</a>	Data Acquired	Apr-18-2013 00:00:00	2013-04-19 15:41:22	\\isaac\labinstruments\Bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	BA2100-4-EQU283	Production	Pending	N/A	-
<input checked="" type="checkbox"/>	9	678760	<a href="#">A29920:CD_1436</a>	<a href="#">133322</a>	Data Acquired	Apr-18-2013 00:00:00	2013-04-19 15:41:22	\\isaac\labinstruments\Bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	BA2100-4-EQU283	Production	Pending	N/A	-
<input checked="" type="checkbox"/>	10	678761	<a href="#">A29921:CD_1437</a>	<a href="#">133323</a>	Data Acquired	Apr-18-2013 00:00:00	2013-04-19 15:41:22	\\isaac\labinstruments\Bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	BA2100-4-EQU283	Production	Pending	N/A	-
<input checked="" type="checkbox"/>	11	678762	<a href="#">A29922:CD_1438</a>	<a href="#">133324</a>	Data Acquired	Apr-18-2013 00:00:00	2013-04-19 15:41:22	\\isaac\labinstruments\Bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	BA2100-4-EQU283	Production	Pending	N/A	-

11 Records

Add Analysis Data





Operation and Maintenance of the Agilent 2100 Bioanalyzer for RNA Samples	
Saved as: LIBPR.0018	Supersedes: Version 8
Version: 9	Page 17 of 19

## Non Controlled Version

\*Note: Controlled Versions of this document are subjected to change without notice

- 1.10. Depending on whether your run is an mRNA assay or a RNA nano assay, different attributes will show up for you to fill out. Once you have finished filling out the information, select “Submit Analysis Data”.

**Enter Bioanalyzer Run Analysis Results:**

RNA Bioanalyzer Analysis Results:

Positions:	1,2,3,4,5,6,7,8,9,10,11
Runs:	133314,133315,133316,133317,133318,133319,133320,133321,133322,133323,133324
Plates:	678752,678753,678754,678755,678756,678757,678758,678759,678760,678761,678762
Libraries:	A29912,A29913,A29914,A29915,A29916,A29917,A29918,A29919,A29920,A29921,A29922
Dilution Factor:	<input style="width: 100%;" type="text" value="5, 5, 5, 5, 5, 4, 4, 4, 4, 4, 5"/>
Agilent Concentration:	<input style="width: 100%;" type="text" value="100, 75, 50, 25, 200, 175, 150, 125, 300, 275, 250"/>
RIN:	<input style="width: 100%;" type="text" value="10, 9, 8, 7, 7, 8, 9, 10, 9, 8, 7"/>

**Submit Analysis Data**



<b>Operation and Maintenance of the Agilent 2100 Bioanalyzer for RNA Samples</b>	
<b>Saved as: LIBPR.0018</b>	<b>Supersedes: Version 8</b>
<b>Version: 9</b>	<b>Page 18 of 19</b>

## Non Controlled Version

\*Note: Controlled Versions of this document are subjected to change without notice

1.11. The results table will look like this. If the volumes are entered for your samples than lims will calculate the amount for you.

Successfully updated Bioanalyzer run analysis data for Run(s) 133314, 133315, 133316, 133317, 133318, 133319, 133320, 133321, 133322, 133323, 133324

Plate Attributes have been added/updated

Experimental Data: Bioanalyzer_Run - RNA																
Run Type	Experiment	Run Status	Validation	Library	Original Source	Plate/Tube	Position	Result File	CHIP ID	Assay Type	Dilution Factor	Agilent Concentration	RIN	Agilent concentration ng uL	RNA Amount ng	
Bioanalyzer_Run	<a href="#">RUN133314</a> <a href="#">A29912-1</a>	Analyzed	Pending	A29912:CD_1408	40856:CD_1408	<a href="#">Pkg78752</a> <a href="#">A29912-1</a>	1	\\isaac\labinstruments\bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	5	100	10	500	10000	<a href="#">Add/Edit Analysis Data</a>
Bioanalyzer_Run	<a href="#">RUN133315</a> <a href="#">A29913-1</a>	Analyzed	Pending	A29913:CD_1429	40857:CD_1429	<a href="#">Pkg78753</a> <a href="#">A29913-1</a>	2	\\isaac\labinstruments\bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	5	75	9	375	7500	<a href="#">Add/Edit Analysis Data</a>
Bioanalyzer_Run	<a href="#">RUN133316</a> <a href="#">A29914-1</a>	Analyzed	Pending	A29914:CD_1430	40858:CD_1430	<a href="#">Pkg78754</a> <a href="#">A29914-1</a>	3	\\isaac\labinstruments\bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	5	50	8	250	5000	<a href="#">Add/Edit Analysis Data</a>
Bioanalyzer_Run	<a href="#">RUN133317</a> <a href="#">A29915-1</a>	Analyzed	Pending	A29915:CD_1431	40859:CD_1431	<a href="#">Pkg78755</a> <a href="#">A29915-1</a>	4	\\isaac\labinstruments\bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	5	25	7	125	2500	<a href="#">Add/Edit Analysis Data</a>
Bioanalyzer_Run	<a href="#">RUN133318</a> <a href="#">A29916-1</a>	Analyzed	Pending	A29916:CD_1432	40860:CD_1432	<a href="#">Pkg78756</a> <a href="#">A29916-1</a>	5	\\isaac\labinstruments\bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	5	200	7	1000	20000	<a href="#">Add/Edit Analysis Data</a>
Bioanalyzer_Run	<a href="#">RUN133319</a> <a href="#">A29917-1</a>	Analyzed	Pending	A29917:CD_1433	40861:CD_1433	<a href="#">Pkg78757</a> <a href="#">A29917-1</a>	6	\\isaac\labinstruments\bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	4	175	8	700	14000	<a href="#">Add/Edit Analysis Data</a>
Bioanalyzer_Run	<a href="#">RUN133320</a> <a href="#">A29918-1</a>	Analyzed	Pending	A29918:CD_1434	40862:CD_1434	<a href="#">Pkg78758</a> <a href="#">A29918-1</a>	7	\\isaac\labinstruments\bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	4	150	9	600	12000	<a href="#">Add/Edit Analysis Data</a>
Bioanalyzer_Run	<a href="#">RUN133321</a> <a href="#">A29919-1</a>	Analyzed	Pending	A29919:CD_1435	40863:CD_1435	<a href="#">Pkg78759</a> <a href="#">A29919-1</a>	8	\\isaac\labinstruments\bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	4	125	10	500	10000	<a href="#">Add/Edit Analysis Data</a>
Bioanalyzer_Run	<a href="#">RUN133322</a> <a href="#">A29920-1</a>	Analyzed	Pending	A29920:CD_1436	40864:CD_1436	<a href="#">Pkg78760</a> <a href="#">A29920-1</a>	9	\\isaac\labinstruments\bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	4	300	9	1200	24000	<a href="#">Add/Edit Analysis Data</a>
Bioanalyzer_Run	<a href="#">RUN133323</a> <a href="#">A29921-1</a>	Analyzed	Pending	A29921:CD_1437	40865:CD_1437	<a href="#">Pkg78761</a> <a href="#">A29921-1</a>	10	\\isaac\labinstruments\bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	4	275	8	1100	22000	<a href="#">Add/Edit Analysis Data</a>
Bioanalyzer_Run	<a href="#">RUN133324</a> <a href="#">A29922-1</a>	Analyzed	Pending	A29922:CD_1438	40866:CD_1438	<a href="#">Pkg78762</a> <a href="#">A29922-1</a>	11	\\isaac\labinstruments\bioanalyzer_Run\BA2100-4\RNA#1550\RNA#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.xad	RNA#1550	Eukaryote Total RNA Nano	5	250	7	1250	25000	<a href="#">Add/Edit Analysis Data</a>

1.12. If you go back to the Bioanalyzer Runs tab, your samples will show up under “Analyzed Runs”.

## Non Controlled Version

\*Note: Controlled Versions of this document are subjected to change without notice

Lab
Bioanalyzer Runs
Admin

---

**Bioanalyzer Runs**

No Waiting Analysis runs found;

---

No Pending runs found;

---

**Analyzed Runs**

Last Week Bioanalyzer Runs - RNA																
Run Type	Experiment	Run Status	Validation	Library	Original Source	Plate/Tube	Position	Result File	CHIP ID	Assay Type	Dilution Factor	Agilent Concentration	RIN	Agilent concentration ng uL	RNA Amount ng	
Bioanalyzer_Run	<a href="#">RUN133314</a> <a href="#">A29912-1</a>	Analyzed	Pending	A29912:CD_1428	40856: CD_1428	<a href="#">Plex78752</a> <a href="#">A29912-1</a>	1	Visaac\labinstruments\Bioanalyzer_Run\BA2100-4RNM#1550RNM#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.iad	RNA#1550	Eukaryote Total RNA Nano	5	100	10	500	10000	<a href="#">Add/Edit Analysis Data</a>
Bioanalyzer_Run	<a href="#">RUN133315</a> <a href="#">A29913-1</a>	Analyzed	Pending	A29913:CD_1429	40857: CD_1429	<a href="#">Plex78753</a> <a href="#">A29913-1</a>	2	Visaac\labinstruments\Bioanalyzer_Run\BA2100-4RNM#1550RNM#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.iad	RNA#1550	Eukaryote Total RNA Nano	5	75	9	375	7500	<a href="#">Add/Edit Analysis Data</a>
Bioanalyzer_Run	<a href="#">RUN133316</a> <a href="#">A29914-1</a>	Analyzed	Pending	A29914:CD_1430	40858: CD_1430	<a href="#">Plex78754</a> <a href="#">A29914-1</a>	3	Visaac\labinstruments\Bioanalyzer_Run\BA2100-4RNM#1550RNM#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.iad	RNA#1550	Eukaryote Total RNA Nano	5	50	8	250	5000	<a href="#">Add/Edit Analysis Data</a>
Bioanalyzer_Run	<a href="#">RUN133317</a> <a href="#">A29915-1</a>	Analyzed	Pending	A29915:CD_1431	40859: CD_1431	<a href="#">Plex78755</a> <a href="#">A29915-1</a>	4	Visaac\labinstruments\Bioanalyzer_Run\BA2100-4RNM#1550RNM#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.iad	RNA#1550	Eukaryote Total RNA Nano	5	25	7	125	2500	<a href="#">Add/Edit Analysis Data</a>
Bioanalyzer_Run	<a href="#">RUN133318</a> <a href="#">A29916-1</a>	Analyzed	Pending	A29916:CD_1432	40860: CD_1432	<a href="#">Plex78756</a> <a href="#">A29916-1</a>	5	Visaac\labinstruments\Bioanalyzer_Run\BA2100-4RNM#1550RNM#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.iad	RNA#1550	Eukaryote Total RNA Nano	5	200	7	1000	20000	<a href="#">Add/Edit Analysis Data</a>
Bioanalyzer_Run	<a href="#">RUN133319</a> <a href="#">A29917-1</a>	Analyzed	Pending	A29917:CD_1433	40861: CD_1433	<a href="#">Plex78757</a> <a href="#">A29917-1</a>	6	Visaac\labinstruments\Bioanalyzer_Run\BA2100-4RNM#1550RNM#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.iad	RNA#1550	Eukaryote Total RNA Nano	4	175	8	700	14000	<a href="#">Add/Edit Analysis Data</a>
Bioanalyzer_Run	<a href="#">RUN133320</a> <a href="#">A29918-1</a>	Analyzed	Pending	A29918:CD_1434	40862: CD_1434	<a href="#">Plex78758</a> <a href="#">A29918-1</a>	7	Visaac\labinstruments\Bioanalyzer_Run\BA2100-4RNM#1550RNM#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.iad	RNA#1550	Eukaryote Total RNA Nano	4	150	9	600	12000	<a href="#">Add/Edit Analysis Data</a>
Bioanalyzer_Run	<a href="#">RUN133321</a> <a href="#">A29919-1</a>	Analyzed	Pending	A29919:CD_1435	40863: CD_1435	<a href="#">Plex78759</a> <a href="#">A29919-1</a>	8	Visaac\labinstruments\Bioanalyzer_Run\BA2100-4RNM#1550RNM#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.iad	RNA#1550	Eukaryote Total RNA Nano	4	125	10	500	10000	<a href="#">Add/Edit Analysis Data</a>
Bioanalyzer_Run	<a href="#">RUN133322</a> <a href="#">A29920-1</a>	Analyzed	Pending	A29920:CD_1436	40864: CD_1436	<a href="#">Plex78760</a> <a href="#">A29920-1</a>	9	Visaac\labinstruments\Bioanalyzer_Run\BA2100-4RNM#1550RNM#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.iad	RNA#1550	Eukaryote Total RNA Nano	4	300	9	1200	24000	<a href="#">Add/Edit Analysis Data</a>
Bioanalyzer_Run	<a href="#">RUN133323</a> <a href="#">A29921-1</a>	Analyzed	Pending	A29921:CD_1437	40865: CD_1437	<a href="#">Plex78761</a> <a href="#">A29921-1</a>	10	Visaac\labinstruments\Bioanalyzer_Run\BA2100-4RNM#1550RNM#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.iad	RNA#1550	Eukaryote Total RNA Nano	4	275	8	1100	22000	<a href="#">Add/Edit Analysis Data</a>
Bioanalyzer_Run	<a href="#">RUN133324</a> <a href="#">A29922-1</a>	Analyzed	Pending	A29922:CD_1438	40866: CD_1438	<a href="#">Plex78762</a> <a href="#">A29922-1</a>	11	Visaac\labinstruments\Bioanalyzer_Run\BA2100-4RNM#1550RNM#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.iad	RNA#1550	Eukaryote Total RNA Nano	5	250	7	1250	25000	<a href="#">Add/Edit Analysis Data</a>

1.13. If you want to edit any of the attributes, you can click on the “Add/Edit Analysis Data” button. For example, if the amount was not calculated for you because the volumes of your samples were not updated, you can edit the volume in LIMS, then go to the Bioanalyzer homepage and select the Run you want to edit under “Analyzed Runs”.

Lab
Bioanalyzer Runs
Admin

---

**Bioanalyzer Runs**

No Waiting Analysis runs found;

---

No Pending runs found;

---

**Analyzed Runs**

Last Week Bioanalyzer Runs - RNA																
Run Type	Experiment	Run Status	Validation	Library	Original Source	Plate/Tube	Position	Result File	CHIP ID	Assay Type	Dilution Factor	Agilent Concentration	RIN	Agilent concentration ng uL	RNA Amount ng	
Bioanalyzer_Run	<a href="#">RUN133314</a> <a href="#">A29912-1</a>	Analyzed	Pending	A29912:CD_1428	40856: CD_1428	<a href="#">Plex78752</a> <a href="#">A29912-1</a>	1	Visaac\labinstruments\Bioanalyzer_Run\BA2100-4RNM#1550RNM#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.iad	RNA#1550	Eukaryote Total RNA Nano	5	100	10	500	10000	<a href="#">Add/Edit Analysis Data</a>
Bioanalyzer_Run	<a href="#">RUN133315</a> <a href="#">A29913-1</a>	Analyzed	Pending	A29913:CD_1429	40857: CD_1429	<a href="#">Plex78753</a> <a href="#">A29913-1</a>	2	Visaac\labinstruments\Bioanalyzer_Run\BA2100-4RNM#1550RNM#1550_Eukaryote Total RNA Nano_DE20901540_2013-04-18_17-20-21.iad	RNA#1550	Eukaryote Total RNA Nano	5	75	9	375	7500	<a href="#">Add/Edit Analysis Data</a>

1.14. You would need to enter the attribute values again and “Submit Analysis Data”.

Operation and Maintenance of the Agilent 2100 Bioanalyzer for RNA Samples	
Saved as: LIBPR.0018	Supersedes: Version 8
Version: 9	Page 20 of 19

## Non Controlled Version

\*Note: Controlled Versions of this document are subjected to change without notice

**Enter Bioanalyzer Run Analysis Results:**

**RNA Bioanalyzer Analysis Results:**

Positions:	1
Runs:	133314
Plates:	678752
Libraries:	A29912
Dilution Factor:	<input style="width: 80%;" type="text" value="10"/>
Agilent Concentration:	<input style="width: 80%;" type="text" value="100"/>
RIN:	<input style="width: 80%;" type="text" value="10"/>

Submit Analysis Data

1.15. The run will be updated with the new values that you have entered and any calculations required will be updated as well.

Successfully updated Bioanalyzer run analysis data for Run(s) 133314

Plate Attributes have been added/updated

Experimental Data: Bioanalyzer\_Run - RNA

Run Type	Experiment	Run Status	Validation	Library	Original Source	Plate/Tube	Position	Result File	CHIP ID	Assay type	Dilution Factor	Agilent Concentration	Agilent concentration RIN	RNA Amount ng uL	RNA Amount ng
Bioanalyzer_Run	<a href="#">RUN133314</a> <a href="#">A29912-1</a>	Analyzed	Pending	A29912_CD_1428	40856: CD_1428	<a href="#">P6670752</a> <a href="#">A29912-1</a>	1	\\visaec\binstruments\Bioanalyzer_Run\BA2100-4R\W#1550\RNA#1550_Eukaryote Total RNA Nano_DE0301540_2013-04-18_17-20-21.vad	RNA#1550	Eukaryote Total RNA Nano	10	100	10 1000	20000	<span style="background-color: #e67e22; color: white; padding: 2px 5px; border: 1px solid black;">Add/Edit Analysis Data</span>