

| | |
|--------------------------------------|-----------------------|
| Plate-based rRNA depletion Version 1 | |
| Document#: LIBPR.0127 | Supersedes: Version 4 |
| Version: 6 | Page 1 of 14 |

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

Plate-based rRNA depletion Version 1

I. Purpose

To remove cytoplasmic rRNAs (nuclear-encoded 5S, 5.8S, 18S and 28S rRNA), and mitochondrial rRNA species (12S and 16S rRNA) from Total RNA using the RNase H protocol from the New England Biolabs (NEB).

II. Scope

All procedures are applicable to the BCGSC Library Core and the Library TechD groups.

III. Policy

This procedure will be controlled under the policies of the Genome Sciences Centre, as outlined in the Genome Sciences Centre High Throughput Production Quality Manual (QM.0001). Do not copy or alter this document. To obtain a copy see a QA associate.

IV. Responsibility

It is the responsibility of all personnel performing this procedure to follow the current protocol. It is the responsibility of the Group Leader to ensure personnel are trained in all aspects of this protocol. It is the responsibility of Quality Assurance Management to audit this procedure for compliance and maintain control of this procedure.

V. References

| Document Title | Document Number |
|---|-----------------|
| NEBNext® rRNA Depletion Kit (Human/Mouse/Rat) | E6310X |

VI. Related Documents

| Document Title | Document Number |
|--|-----------------|
| Operation and Maintenance of the Agilent 2100 Bioanalyzer for DNA samples | LIBPR.0017 |
| Quantifying DNA samples using the Qubit Fluorometer | LIBPR.0030 |
| 96-well DNA Quantification using the dsDNA Quant-iT High Sensitivity Assay Kit and VICTOR ³ V | LIBPR.0108 |
| Operation and Maintenance of the Agilent 2100 Bioanalyzer for RNA Samples | LIBPR.0018 |
| Operation and Maintenance of the LabChipGX for RNA samples | LIBPR.0052 |

| | |
|--------------------------------------|-----------------------|
| Plate-based rRNA depletion Version 1 | |
| Document#: LIBPR.0127 | Supersedes: Version 4 |
| Version: 6 | Page 2 of 14 |

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

| Document Title | Document Number |
|--|-----------------|
| using the HT RNA Assay Total RNA Normalization on the Hamilton Nimbus | LIBPR.0121 |

VII. Safety

All Laboratory Safety procedures will be complied with during this procedure. The required personal protective equipment includes a laboratory coat and gloves. See the material safety data sheets (MSDS) for additional information.

VIII. Materials and Equipment

| Name | Supplier | Number |
|--|-------------------|------------|
| Fisherbrand Textured Nitrile gloves - large | Fisher | 270-058-53 |
| RNAse Zap | Ambion | 9780 |
| Ice bucket – Green | Fisher | 11-676-36 |
| Wet ice | In house | N/A |
| RNAse free 1.5 mL eppendorf tube | Ambion | 12400 |
| Gilson P2 pipetman | Mandel | GF-44801 |
| Gilson P10 pipetman | Mandel | GF-44802 |
| Gilson P20 pipetman | Mandel | GF23600 |
| Gilson P200 pipetman | Mandel | GF-23601 |
| Gilson P1000 pipetman | Mandel | GF-23602 |
| Mandel P200 DF200 tips | Mandel | GF-F171503 |
| Mandel P1000 DF1000 tips | Mandel | GF-F171703 |
| VX-100 Vortex Mixer | Rose Scientific | S-0100 |
| 200µL Rainin tips | Rainin | RT-L200F |
| 20µL Rainin tips | Rainin | RT-L10F |
| 200µL Pipet-Lite | Rainin | L12-200 |
| 20µL Pipet-Lite | Rainin | L12-20 |
| 1250 µL pipette tip, 96tips/rack, 480 tips/cs filter sterile | Mandel Scientific | TM-4445 |
| Pipette-VIAFLO 8 Channel | Mandel Scientific | TM-4124 |
| Large Kimwipes | Fisher | 06-666-117 |
| Black ink permanent marker pen | VWR | 52877-310 |
| Bench Coat (Bench Protection Paper) | Fisher | 12-007-186 |
| Small Autoclave waste bags 10”X15” | Fisher | 01-826-4 |
| DNaseI Amplification Grade 100U | Invitrogen | 18068-015 |
| DEPC water | Ambion | 9922 |
| Mini-centrifuge | Eppendorf | 5417R |

| | |
|--------------------------------------|-----------------------|
| Plate-based rRNA depletion Version 1 | |
| Document#: LIBPR.0127 | Supersedes: Version 4 |
| Version: 6 | Page 3 of 14 |

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

| | | |
|--|----------------------|--------------|
| Thermo Scientific 0.2mL Ultra Rigid Skirted 96-well | Thermoscientific | AB1000-150s |
| Deep-well, 96-well, 1.2 mL, U bottom, low pro, 50/cs | Fisher Scientific | AB1127 |
| NEBNext® rRNA Depletion Kit (Human/Mouse/Rat) | NEB | E6310X |
| RNA MagClean DX | ALINE Biosciences | C-1005 |
| Maxima H Minus First Strand cDNA Synthesis Kit | Thermo-Fisher | K165B001 |
| Second Strand cDNA Synthesis Kit | Invitrogen | A26943 |
| Actinomycin D 10mg/mL | In House | |
| Actinomycin D 5mg | MJS Biolynx | ENZGR3000005 |
| dNTP Blend, 12.5mM with dUTP, 1mL | GeneAmp | N8080270 |
| PCR Clean DX (ALINE Beads) | ALINE Biosciences | C-1003-450 |
| Sterile Filtered Conductive 50µL Tips in Frames | Hamilton | 235979 |
| Sterile Filtered Conductive 300µL Tips in Frames | Hamilton | 235938 |
| Tape Pads | Qiagen | 19570 |
| Foil Tape | VWR | 60941-126 |
| MJ Research Tetrad PTC-225 Thermal Cycler | MJ Research | 8252-30-1004 |
| G127187 DNA Engine Tetrad 2 Peltier Thermal Cycler | MJ Research | 10177BD |
| 2100 Electrophoresis Bioanalyzer Instrument | Agilent | G2939AA |
| LabChip GXII | Perkin Elmer/Caliper | 124582 |
| NIMBUS Liquid Handling Workstation | Hamilton | |


IX. GENERAL GUIDELINES

1. General guidelines and input material

The following is a flow chart depicting the various steps in this protocol that will be described in detail below:

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*



1.1 The recommended input material for this procedure is

| Input material | Input amount |
|-------------------------|-----------------|
| FFPE Total Nucleic Acid | 200-400ng DNA |
| FFPE Total RNA | maximum 2µg RNA |

| | |
|--------------------------------------|-----------------------|
| Plate-based rRNA depletion Version 1 | |
| Document#: LIBPR.0127 | Supersedes: Version 4 |
| Version: 6 | Page 5 of 14 |

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

| | |
|------------|-----------|
| Intact RNA | 120ng RNA |
|------------|-----------|

For FFPE RNA, based on smear analysis on Agilent / Caliper, 100-5000nt RNA fragments should constitute >70% of the total RNA fragments. Consult with your supervisor if the % total is <70% for this range. Input volume to be requested from collaborators is a maximum of 20µL to have suitable concentration for QC. The actual input volume for the first reaction is 35µL/well in DEPC H2O in a 96-well plate. This protocol was validated at the GSC with mouse and human RNA.

- 1.2 The positive control for the ribodepletion procedure is 120ng Universal Human Reference RNA (UHR or FG031) and the negative control is DEPC H2O. An additional control that should be added when cDNA synthesis is being set up is 12ng UHR from the same immediate stock as the one that went through the ribodepletion was derived from. Make a dilution of the original stock so that you have 137ng in 40µL. Of this immediate stock, take 35µL for the ribodepletion control. Take 3.5µL for the non-ribodepleted control when you set up the cDNA reaction (to be topped up with DEPC H2O). This would serve as a non-ribodepleted control. The ribodepleted UHR should give a cDNA yield of <100% of the cDNA level from non-ribodepleted UHR. Check with your supervisor if the yield is >100%. The calculations are based on agilent smear analysis (mass) within 50b-9kb range. The non-ribodepleted control will also serve as a reliable positive control for cDNA synthesis and to ensure that shearing worked later in the library construction protocol. The ribodepleted UHR control and the samples are barely detectable before iPCR.
- 1.3 The beads (**RNA MagClean DX**) to be used in this protocol are different from those used for other protocols as they are certified to be RNase-free.
- 1.4 Ensure proper personal protective equipment is used when handling sample plates, reagents and equipment. Treat everything with, and adhere to, strict RNA handling techniques.
- 1.5 Wipe down the assigned workstation, pipetman, tip boxes and small equipment with RNAse Zap (Ambion) followed by DEPC-treated water. Ensure you have a clean working surface before you start.
- 1.6 Double check the QA release and/or expiry date of each reagent and enzyme.
- 1.7 Reactions in plates should never be vortexed and plate covers are never to be re-used.
- 1.8 Retrieve and thaw all reagents at room temperature. Once thawed, pulse-vortex, quick spin and keep reagents on ice. Enzymes should be left in the freezer until ready to use.
- 1.9 Ensure the waste bag for the Nimbus is empty.

| | |
|--------------------------------------|-----------------------|
| Plate-based rRNA depletion Version 1 | |
| Document#: LIBPR.0127 | Supersedes: Version 4 |
| Version: 6 | Page 6 of 14 |

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

- 1.10 Brews are prepared and dispensed manually. Note that what is dispensed into the final brew plate is the actual volume of the specified volume for each step below without any dead volume.
- 1.11 Beads are manually dispensed into wells in a 1.2 ml plate (AB1127). 20µL dead volume should be factored in, which means you need to dispense 110µL in the final bead source plate that will be placed on the Nimbus. Ensure that you include 25ml dead volume for 70% ethanol and DEPC H2O on top of what is required for the actual washes and elution, respectively.
- 1.12 The Nimbus adds sample or sample in a given reaction to the brew and bead cleanups are performed on Nimbus. Follow the prompts and lay out from the Nimbus programs to execute a particular step.
- 1.13 The Nimbus mixes at 80% of total volume 10 times.
- 1.14 Note that where it is specified that you proceed immediately to the next step, plates can be briefly placed on ice (not more than 30min) in the case of emergency. Make an active attempt to proceed as immediate as possible.

X. PROCEDURE

Note: All version numbers for Nimbus protocols have been removed on this document. They are present when running the protocol. If you are unsure which version to use, consult your supervisor.

1. Upstream Preparation

- 1.1 In the case that LIMS is down, brew calculators can be located on the worksheet listed below. Enter the number of samples to be processed and print the Ribodepletion Worksheet located in:

| |
|---|
| R:\Library Core\Work Sheets and Calculators\Ribodepletion |
|---|

- 1.2 Retrieve the plate containing Total RNA. If stored in -80°C, thaw it on ice followed by a quick spin at 4°C, 700g for 1min. Place it on ice.

2. Input QC:

- 2.1 For FFPE total nucleic acid, it is recommended that contaminating gDNA is first quantified using Qubit (LIBPR.0030) or Quant-IT (LIBPR.0108). RNA QC can be skipped at this stage (200-400ng gDNA should generally give more than 200ng RNA after DNase treatment). The RNA/DNA mixture entering this DNase treatment should contain <400ng gDNA. Normalization to 200-400ng gDNA should be performed using Nimbus (LIBPR.0121) or manually if deemed appropriate.

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

- 2.2 For FFPE RNA or intact RNA extracted using protocols such as Qiagen's where RNA is separated from gDNA, RNA is quantified using Agilent RNA Nano (LIBPR.0018) or Caliper RNA (LIBPR.0052) assays. Normalize intact RNA samples to 120ng and FFPE RNA to maximum 2µg using Nimbus (LIBPR.0121) or manually if deemed appropriate.

| |
|---|
| LIBPR.0121 Total RNA Normalization on the Hamilton Nimbus |
|---|

3. 1st DNase I Treatment: Remove contaminating gDNA in samples

All samples should be DNase treated as follows (regardless of whether or not they were treated by collaborators). Samples should be in 35µL.

- 3.1 For intact RNA extracted using protocols such as Qiagen's where RNA is separated from gDNA, the reaction set up for 1 reaction is shown below. Use LIMS to generate the brew calculator for this step. If LIMS is not working, use the excel worksheet in the path above.

| Solution | µL (per 1rxn) |
|------------------------------|---------------|
| RNA | 35 |
| DEPC H2O | 9 |
| 10X DNase Buffer | 5 |
| DNase I Enzyme | 1 |
| Total Reaction volume | 50 |

| |
|--------------------------------------|
| 1 st DNase Mix (15 µL) |
|--------------------------------------|

For FFPE RNA or FFPE Total Nucleic acid, the reaction set up for 1 reaction is shown below:

| Solution | µL (per 1rxn) |
|------------------|---------------|
| RNA | 35 |
| DEPC H2O | 5 |
| 10X DNase Buffer | 5 |
| DNase I Enzyme | 5 |

| |
|--------------------------------------|
| 1 st DNase Mix (15 µL) |
|--------------------------------------|

| | |
|--------------------------------------|-----------------------|
| Plate-based rRNA depletion Version 1 | |
| Document#: LIBPR.0127 | Supersedes: Version 4 |
| Version: 6 | Page 8 of 14 |

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

| | |
|-----------------------|----|
| Total Reaction volume | 50 |
|-----------------------|----|

- 3.2 Make sure to mix each reagent well and quick spin before adding to the brew. Add enzyme last.
- 3.3 Prepare the reaction brew in a non-stick tube and check off reagents as they are added on the worksheet. Mix the brew by repeated pulse-vortexing followed by a quick spin.

LIMS Calculator (intact RNA): RBD_1st DNase_1U

LIMS Calculator (FFPE RNA or Total Nucleic acid): RBD_1st DNase_5U

- 3.4 Log into Nimbus Program as follows to add the RNA to 15µL of the DNase mix:

Open file: Production > toggle to workflow > Ribodepletion > Ribodepletion Scheduler.wfl > Brew Additions > 1st DNase

- 3.5 After Nimbus program completion, seal the plates and quick spin at 4°C for 1 minute. Inspect the reaction plates for any variations in volume.
- 3.6 Incubate for 15 min at room temperature. *Note: This is NOT a safe stopping point; proceed to the next step immediately.*
- 3.7 Log into Nimbus Program as follows to add the RNA to 5µL of 25mM EDTA:

Open file: Production > toggle to workflow > Ribodepletion > Ribodepletion Scheduler.wfl > Brew Additions > EDTA

- 3.8 After Nimbus program completion, seal the plates and quick spin at 4°C for 1 minute. Inspect the reaction plates for any variations in volume.
- 3.9 In the Tetrad thermocycler, incubate the plate at 65°C for 10 min. During incubation, dispense 110 µL of Beads per well (1.2ml plate) in preparation for the subsequent cleanup.

TETRAD: RBD > RBD_65

- 3.10 After incubation, remove the plate and spin it at 4°C, 700g for 1min. *Note: This is NOT a safe stopping point; proceed to the next step immediately.*
- 4. Post-1st DNase I treatment Bead clean up (use RNA MagClean DX beads)**

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

- 4.1 The input volume for this step is 50µL per well. Note that even though the previous reaction is in 55µL, you will be prompted to transfer 50µL. The supernatant after 1st clearing from the bead cleanup may or may not be used as miRNA fraction. The reason for transferring 50µL is to ensure consistent size selection in the case that miRNA fraction is required. Check with your supervisor if the miRNA fraction needs to be processed.
- 4.2 The Nimbus will perform the cleanup of the 1st DNase reaction using beads as follows:

| Reaction | *Bead Vol manually dispensed (µL) | Bead Vol added by Nimbus (µL) | Bead Binding Time (mins) | Magnet Clearing Time (mins) | 2X 70% EtOH* Wash Vol (µL) | Magnet Airdry Time (mins) | DEPC Elution Vol (µL) | Elution time (mins) | Magnet Elution time (mins) | Transfer Vol (µL) |
|--------------------|-----------------------------------|-------------------------------|--------------------------|-----------------------------|----------------------------|---------------------------|-----------------------|---------------------|----------------------------|-------------------|
| 1st DNase Reaction | 110 | 90 | 15 | 7 | 180 | 5 | 15 | 3 | 2 | 15 |

*Must be at Room Temp for a minimum of 30mins before usage; failure to do so would result in a decrease in yield

- 4.3 Log into Nimbus Program as follows:

Open file: Production > toggle to workflow > Ribodepletion > Ribodepletion Scheduler.wfl > Bead Clean > RNA Bead Cleanup > RBD Post 1st DNase

If the miRNA fraction is to be kept, click on the box when asked if you would like to keep the supernatant. Follow the prompts to place a 1.2mL Deep-well plate in the destination position. Once supernatant has been transferred, remove plate, seal with foil tape, label as ‘supernatant’, and store at -80°C. Continue with bead clean as instructed by the Nimbus.

This is a safe stopping point. If needed, the plate containing the bead cleaned RNA can be stored at -80°.

- 4.4 For FFPE Total Nucleic Acid and FFPE RNA, DNase-treated RNA should be quantified using Agilent RNA Nano (LIBPR.0018) or Caliper RNA (LIBPR.0052). Skip this QC for intact RNA. Based on smear analysis, 100-5000nt RNA fragments should constitute >70% of the total RNA fragments. Consult with your supervisor if the % total is <70% for this range. Confirm with your supervisor on which assay to use.
- 4.4.1 For Caliper QC on Standard Assay, use 2µL of sample and 46µL of made up HT RNA Caliper buffer. The Quantitative range for the Standard Assay is 25-250ng/µL.

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

4.4.2 For Caliper QC on HiSens Assay, use 2µL of sample diluted with 4µL DEPC water (6µL total) and 19µL of made up HT RNA Caliper buffer. Also, select “AB1000_2.5mm” as the total volume is lower for the Caliper plate. The Quantitative range for the HiSens Assay is 5-50ng/µL.

4.5 For Caliper QC setup on the Nimbus:

Open file: Production > toggle to workflow > Ribodepletion > Ribodepletion Scheduler.wfl > Caliper Transfer > Caliper Transfer > Standard or High Sensitivity

4.6 DNase-treated RNA from FFPE Total nucleic Acid and FFPE RNA should be normalized to 120ng at this stage using Nimbus according to LIBPR.0121.

5. Hybridization: Annealing of rRNA probes to target rRNAs

Note: all reagents for steps 5, 6, and 7 including the DNase I in Step 7 are from the NEBNext rRNA depletion kit.

5.1 The reaction set up for 1 reaction is as follows:

| Solution | 1 rxn (µL) |
|-----------------------------|------------|
| DNase-treated RNA | 12 |
| Depletion solution (probes) | 1 |
| Hybridization Buffer | 2 |
| Total volume | 15 |

} Hyb Mix (3 µL)

5.2 Make sure to mix each reagent well and quick spin before adding to the brew.

5.3 Prepare the brew and check off reagents as they are added. Mix the brew by repeated pulse-vortexing followed by a quick spin.

LIMS Calculator: RBD_Hybridization

5.4 Log into Nimbus Program as follows to add the RNA to 3 µL of the Hyb mix:

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

Open file: Production > toggle to workflow > Ribodepletion > Ribodepletion Scheduler.wfl > Brew Additions > Hybridization

5.5 After Nimbus program completion, seal the plates and quick spin at 4°C for 1 minute. Inspect the reaction plates for any variations in volume.

5.6 In the MJ Research Gradient Tetrads thermocycler, incubate the plate as follows:

TETRAD: RNAHYB

- 95°C 2 min
- 95°C, -0.1°C at 1 sec down to 22°C (730 cycles)
- 22°C 5min

Attend this immediately after the incubation is done (within 30min).

5.7 After thermo-cycler program is finished, remove the plate and spin it at 4°C, 700g for 1min. *Note: This is NOT a safe stopping point; proceed to the next step immediately.*

6. RNase H digestion: Degradation of RNA of in rRNA/ DNA probe hybrid

6.1 The reaction set up for 1 reaction is as follows:

| Solution | 1 rxn (µL) |
|-------------------------|------------|
| RNA in hyb reaction | 15 |
| Nuclease free water | 1 |
| RNase H Reaction Buffer | 2 |
| RNase H | 2 |
| Total volume | 20 |

} RNase H mix (5 µL)

6.2 Make sure to mix each reagent well and quick spin before adding to the brew.

6.3 Prepare the brew and check off reagents as they are added. Mix the brew by repeated pulse-vortexing followed by a quick spin.

LIMS Calculator: RBD_ RNaseH

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

6.4 Log into Nimbus Program as follows to add the RNA to 5 µL of the RNase H mix:

Open file: Production > toggle to workflow > Ribodepletion > Ribodepletion Scheduler.wfl > Brew Additions > RNaseH

6.5 After Nimbus program completion, seal the plates and quick spin at 4°C for 1 minute. Inspect the reaction plates for any variations in volume.

6.6 In the Tetrad thermocycler, incubate the plate at 37°C for 30 min (Heated Lid 47°C).

TETRAD: RBD>RBD_37

6.7 After incubation, remove the plate and spin it at 4°C, 700g for 1min. *Note: This is NOT a safe stopping point; proceed to the next step immediately.*

7. 2nd DNase I treatment: Degradation of rRNA DNA probes

7.1 The reaction set up for 1 reaction is as follows (all reagents for this reaction should be from the Ribodepletion kit E6310X).

| Solution | 1 rxn (µL) |
|----------------------------------|------------|
| RNA in RNase H reaction | 20 |
| Nuclease free water | 22.5 |
| DNase Buffer (from NEBNext Kit) | 5 |
| DNase I (from NEBNext Kit) | 2.5 |
| Total volume | 50 |

DNase mix (30 µL)

7.2 Make sure to mix each reagent well and quick spin before adding to the brew.

7.3 Prepare the brew and check off reagents as they are added. Mix the brew by repeated pulse-vortexing followed by a quick spin.

LIMS Calculator: RBD_ 2nd DNase

7.4 Log into Nimbus Program as follows to add the RNA to 30 µL of the DNase mix:

Open file: Production > toggle to workflow > Ribodepletion >

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

Ribodepletion Scheduler.wfl > Brew Additions > 2nd DNase

- 7.5 After Nimbus program completion, seal the plates and quick spin at 4°C for 1 minute. Inspect the reaction plates for any variations in volume.
- 7.6 In the Tetrad thermocycler, incubate the plate at 37°C for 30 min (Heated Lid 47°C). During incubation dispense 110 µL of Beads per well (1.2ml plate) in preparation for the subsequent cleanup.

TETRAD: RBD>RBD_37

- 7.7 After incubation, remove the plate and spin it at 4°C, 700g for 1min. *Note: This is NOT a safe stopping point; proceed to the next step immediately.*
8. **Post-2nd DNase treatment Bead clean up (use RNA MagClean DX beads)**

- 8.1 The input volume for this step is 50µL per well.
- 8.2 The Nimbus will perform the cleanup of the 2nd DNase reaction using beads as follows:

| Reaction | *Bead Vol manually dispensed (µL) | Bead Vol added by Nimbus (µL) | Bead Binding Time (mins) | Magnet Clearing Time (mins) | 2X 70% EtOH* Wash Vol (µL) | Magnet Airdry Time (mins) | DEPC Elution Vol (µL) | Elution time (mins) | Magnet Elution time (mins) | Transfer Vol (µL) |
|--------------------------------------|-----------------------------------|-------------------------------|--------------------------|-----------------------------|----------------------------|---------------------------|-----------------------|---------------------|----------------------------|-------------------|
| 2nd DNase Reaction | 110 | 90 | 15 | 7 | 180 | 5 | 37 | 3 | 2 | 37 |

*Must be at Room Temp for a minimum of 30mins before usage; failure to do so would result in a decrease in yield

- 8.3 Log into Nimbus Program as follows:

Open file: Production > toggle to workflow > Ribodepletion > Ribodepletion Scheduler.wfl > Bead Cleans > RNA Bead Cleanup > RBD Post 2nd DNase

This is a safe stopping point. If needed, the plate containing the bead cleaned RNA can be stored at -80°C. Otherwise, proceed to cDNA synthesis as described in 96-well Plate-based Strand-specific cDNA Synthesis using Maxima H Minus on Hamilton NIMBUS-LIBPR.0132.

Non-Controlled Version

**Note: Controlled Versions of this document are subjected to change without notice*

Appendix A: LIMS SOP

1. Start of Plate Library Construction
2. Bioanalyzer Run / Caliper Run – if working with FFPE RNA or FFPE total nucleic acid, QC samples on Agilent after 1st DNase treatment. Create Bioanalyzer Run – QC Category: Total RNA QC
3. Ribodepletion-1st DNase Treatment: pass barcode to supervisor
4. Ribodepletion: Use the barcode your supervisor passes to you.
5. Pass barcode to supervisor to add control for cDNA synthesis.